

Микроскопы сверхвысокого разрешения Lattice SIM 3, Lattice SIM 3, Elyra 7 with Lattice SIM

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: zsf@nt-rt.ru || сайт: <https://zeiss.nt-rt.ru/>

Super-resolution imaging across scales



ZEISS Lattice SIM Family

Full Access to Super-Resolution Imaging
for all Research Areas




Seeing beyond

ZEISS Lattice SIM Family


Full Access to Super-Resolution Imaging for all Research Areas

Using microscopy to visualize biological structures provides insights into function. When imaging fixed structures, acquisition settings can be optimized for spatial resolution. However, when capturing dynamic events in living samples, higher acquisition speeds and low-light conditions must be balanced with resolution. The ZEISS Lattice SIM family balances sample size, imaging speed, and super-resolution capabilities based on your application – from outstanding optical sectioning of tissues and developing organisms to high-speed imaging of living cells to resolution excellence at the molecular level.



ZEISS Lattice SIM 3


Your fast optical sectioning solution for studying developing organisms and tissue microstructures



Spheroid stained for mitochondria (MitoTracker Green) and nuclei (NucRed Live 647)

ZEISS Lattice SIM 5


Your live imaging system for uniform super-resolution in all spatial dimensions



Cos-7 cells stained for microtubules (anti-tubulin Alexa Fluor 488, cyan) and actin (Phalloidin Alexa Fluor 561, orange)

ZEISS Elyra 7 with Lattice SIM

Your complete super-resolution system with unprecedented resolution down to molecular details




Dual-color 2D STORM of Cos-7 cells stained for microtubules (anti-tubulin-Alexa Fluor 647, magenta) and mitochondria (anti-TOMM20-CF568, yellow)

ZEISS Lattice SIM 3

Reveal Cellular Behavior and Inter-Cellular Dynamics

ZEISS Lattice SIM 3 is designed to meet the requirements of multicellular samples, such as: developing organisms, organoids, 3D cell cultures, and tissue sections. Optimized for use with objectives from 10x to 40x, ZEISS Lattice SIM 3 exploits the full potential of the SIM Apotome technology: fast optical sectioning at superior quality, large fields of view with access to smaller regions of interest, near-isotropic resolution, and the gentlest super-resolution imaging possible. Additionally, Lattice SIM imaging and SIM² image reconstruction give you super-resolution imaging down to 140 nm.



Drosophila embryo stained for Fasciclin II (color-coded depth projection) and HRP (cyan) labelling the nervous system, imaged in SIM Apotome mode, Objective: Plan-Apochromat 20x/0.8 air. Image courtesy of Ines Hahn, University of York, UK.

Capture entire model organisms and tissue sections

ZEISS Lattice SIM 3 fully leverages SIM Apotome technology, to provide the most outstanding optical sectioning at large fields of view with near-isotropic resolution. ZEISS Lattice SIM 3 is your system of choice for fast imaging of larger volumes, such as 3D model organisms, embryos, organoids, or tissue sections. Whether you work with living or fixed samples, ZEISS Lattice SIM 3 provides access to structured illumination microscopy of multicellular organisms with superior penetration depth.

Acquire super-resolution images as fast and gentle as widefield images

Choose between the standard SIM Apotome imaging mode for the highest available resolution (5 phase images required) or the imaging mode with reduced phases for slightly lower resolution but significantly increased speed and gentleness (only 3 phase images required). Combine SIM Apotome with the Leap mode to significantly speed up super-resolution acquisition. SIM Apotome makes even lossless acquisition possible, meaning for every reconstructed image just one raw image is needed.



Go from a large-field overview to the super-resolution details

For large sample experiments, ZEISS Lattice SIM 3 offers the most advantageous combination of a large field of view and super-resolution imaging. SIM Apotome mode in combination with SIM² image reconstruction enables lateral super-resolution down to 140 nm with superior optical sectioning and sensitivity. Additionally, imaging in Lattice SIM mode with a ZEISS 25x multi-immersion objective and subsequent SIM² processing provides similar lateral resolutions with larger fields of view and more flexible adaptation to the refractive index of your sample.

ZEISS Lattice SIM 5

Reveal the Vibrant Sub-Organelle Network of Life

ZEISS Lattice SIM 5 has been optimized for single cell imaging as well as capturing subcellular structures and their dynamics. Powered by the Lattice SIM technology and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 provides you with outstanding super-resolution capabilities down to 60 nm in both living and fixed cells. Additionally, you can choose SIM Apotome imaging mode and a low-magnification objective to achieve fast overview images of your sample before zooming into super-resolution details.



Murine brain expressing the neuronal marker Thy1-eGFP was imaged in SIM Apotome and Lattice SIM modes over a Z stack range of 170 µm. Objective for overview image (left): Plan-Neofluar 10x. The ZEN Connect project combines data sets recorded with 10x SIM Apotome, 25x SIM Apotome, 40x SIM Apotome and 63x Lattice SIM. The volume rendering on the right-hand side shows a subset of the 63x Lattice SIM data set. Objective: Plan-Apochromat 63x/1.4 Oil. Sample courtesy of Herms Lab (MCN, University of Munich, Germany).

Capture highly dynamic processes

Equipped with Lattice SIM illumination and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 raises structured illumination microscopy to a new level. You will always achieve the best possible results, even when using lower light exposures to protect living specimens. Double the conventional SIM resolution and discriminate the finest subcellular structures that are no more than 60 nm apart. The light-efficient Lattice SIM technology provides the gentlest imaging of living and fixed specimens, giving you high temporal resolution with up to 255 fps.

Optimize to the needs of living samples

The flexibility of ZEISS Lattice SIM 5 allows you to balance the needs of your experiment by prioritizing resolution, speed, or by finding the right balance in between. Use the photon budget to enhance lateral resolution well below 100 nm or reduce the number of required raw images to boost acquisition speed and gentleness. ZEISS Lattice SIM 5 has a number of options for reducing raw images which allows you to select for the best acquisition settings that target your desired spatial and temporal resolution.


Get more reliable experiment results

ZEISS Lattice SIM 5 with SIM² comes with outstanding out-of-focus light suppression, giving you the sharpest sectioning in widefield microscopy, even for highly scattering samples. The SIM² image reconstruction uses a special SIM point spread function to robustly reconstruct all structured-illumination-based acquisition data of your ZEISS Lattice SIM 5 with minimal image artifacts – for both living and fixed samples. Rest assured knowing that you are basing your experimental conclusions on reproducible data generated from a powerful and proven algorithm.

ZEISS Elyra 7 with Lattice SIM

Reveal Life Across Scales – Down to Molecular Details

ZEISS Elyra 7 includes a wealth of microscopy techniques to meet your experimental needs across scales, optimally matching resolution, speed, and sensitivity requirements to your demanding samples. Employ SIM Apotome for fast optical sectioning, Lattice SIM for super-resolution imaging, SIM² image reconstruction for resolution excellence down to 60 nm, as well as SMLM and TIRF for investigations at the molecular level. You can combine these techniques to multiply the insights from your specimen and to correlate the acquired data.



Single-color 3D STORM of *Xenopus laevis* A6 cells stained for nuclear pore marker GP210 (Georg Krohne, University of Würzburg) and F(ab)2 Alexa Fluor 647 secondary antibody. Objective: alpha Plan-Apochromat 63x/ 1.46 oil.

Observe life's finest details

ZEISS Elyra 7 was designed for the understanding of biological processes that require the observation of living samples in low light and high spatiotemporal resolution. Observe cellular, sub-cellular, and sub-organelle structures in living specimens in 2D and 3D over time. Whether you are interested in cytoskeletal dynamics, mitochondrial fusion and fission, or budding of the endoplasmic reticulum, it gives you both live cell compatibility and super-resolution.

3D imaging at molecular resolution

With Elyra 7 you can use single-molecule localization microscopy (SMLM) techniques such as PALM, dSTORM and PAINT to achieve lateral localization precision of 10 – 20 nm. The ZEN software will seamlessly perform the image reconstruction of your data. In addition, Elyra 7 provides you with 3D SMLM mode based on PRILM technology. The PSF is reshaped for encoding the Z position so while acquiring only one plane, you get volume information of 1.4 μm depth at 20 – 40 nm axial resolution. Thus, you can acquire 3D data with consistent molecular precision.

Correlative microscopy within the same system

Correlative microscopy, where the same area of the sample is imaged using different techniques, has become an important tool for biological research. Imaging with Elyra 7 can be combined seamlessly with LSM 980 Airyscan or ZEISS electron microscopy solutions. Nevertheless, Elyra 7 itself contains three different imaging modalities – SIM Apotome, Lattice SIM and SMLM – providing the possibility to combine them when needed for sample scales varying by orders of magnitude.

ZEISS Lattice SIM Family

Balance Sample Size, Imaging Speed, and Super-Resolution Capabilities Based on the Needs of Your Application



	Lattice SIM 3	Lattice SIM 5	Elyra 7 with Lattice SIM
Target applications	<ul style="list-style-type: none">■ 3D cell culture models■ Organoid development■ Embryonic development■ Tissue sections	<ul style="list-style-type: none">■ Fast dynamics in living cells■ Organelle and protein interaction■ Structural super-resolution imaging of fixed samples	<ul style="list-style-type: none">■ Ultrastructure of organelles and molecular assemblies■ Single molecule localization microscopy (SMLM)■ Molecule interaction■ Protein localization
SIM Apotome module	Two different grating frequencies for best match to objective and wavelength	One grating	One grating
Lattice SIM	One grating	Five different grating frequencies for optimal matching of illumination pattern to laser wavelength and objective lens	
Multi-tracking	Freely configurable change of gratings (SIM Apotome); One common grating (Lattice SIM)	Freely configurable change of gratings (Lattice SIM); One common grating (SIM Apotome)	
Leap mode	✓	✓	✓
Burst mode	✓	✓	✓
TIRF illumination			✓
3D-PALM			✓
Resolution			
SIM/SIM² Apotome	XY: down to 320 / 265 nm for 25×		
Lattice SIM/SIM²	XY: down to 210 nm / 140 nm with LD LCI Plan-Apochromat 25× / 0.8 Imm Corr DIC	XY: down to 120 nm / 60 nm with Plan-Apochromat 63× / 1.40 Oil DIC	
Objective lenses			
SIM Apotome	Plan-Apochromat 40×/1.4 Oil; C-Apochromat 40×/1.2 W; LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC; Plan-Apochromat 20×/0.8 Air; EC Plan-Neofluar 10×/0.3 Air		
Lattice SIM	LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC	Plan-Apochromat 63×/1.40 Oil DIC; C-Apochromat 63×/1.20 W Corr; alpha Plan-Apochromat 63×/1.46 Oil	Plan-Apochromat 63×/1.40 Oil DIC; C-Apochromat 63×/1.20 W Corr; alpha Plan-Apochromat 63×/1.46 Oil; alpha Plan-Apochromat 100×/1.57 Oil-HI DIC Corr
SMLM			alpha Plan-Apochromat 63×/1.46 Oil (2D/3D PALM); Plan-Apochromat 63×/1.4 Oil DIC (3D PALM); C-Apochromat 63×/1.2 W Corr DIC (3D-PALM); Objectives with NA >= 1.46 suitable for TIRF and HILO illumination
TIRF			alpha Plan-Apochromat 100×/1.46 Oil DIC; alpha Plan-Apochromat 100×/1.57 Oil-HI DIC Corr (2D-PALM);
Camera	Up to two CMOS cameras (ZEISS AxioCam 820 mono)	Up to two CMOS cameras (ZEISS AxioCam 820 mono) or up to two sCMOS cameras (Hamamatsu ORCA-Fusion BT)	Up to two sCMOS cameras (Hamamatsu ORCA-Fusion BT)
Multi-color imaging	Detection of up to four different fluorescent labels (sequential detection); simultaneous dual-color detection with Duolink		

Revealing cellular behavior and inter-cellular dynamics



ZEISS Lattice SIM 3

Your Fast Optical Sectioning Solution for Studying Developing Organisms and Tissue Microstructures



Seeing beyond

Your Fast Optical Sectioning Solution for Studying Developing Organisms and Tissue Microstructures

› In Brief

› The Advantages

› The Applications

› The System

› Technology and Details

› Service


The ZEISS Lattice SIM family

Using microscopy to visualize biological structures provides insights into function. When imaging fixed structures, acquisition settings can be optimized for spatial resolution. However, when capturing dynamic events in living samples, higher acquisition speeds and low-light conditions must be balanced with resolution. The ZEISS Lattice SIM family balances sample size, imaging speed, and super-resolution capabilities based on your application – from outstanding optical sectioning of tissues and developing organisms to high-speed imaging of living cells to resolution excellence at the molecular level.

ZEISS Lattice SIM 3

ZEISS Lattice SIM 3 is designed to meet the requirements of multicellular samples, such as: developing organisms, organoids, 3D cell cultures, and tissue sections. Optimized for use with objectives from 10× to 40×, ZEISS Lattice SIM 3 exploits the full potential of the SIM Apotome technology: fast optical sectioning at superior quality, large fields of view with access to smaller regions of interest, near-isotropic resolution, and the gentlest super-resolution imaging possible. Additionally, Lattice SIM imaging and SIM² image reconstruction give you super-resolution imaging down to 140 nm.

With ZEISS Lattice SIM 3, not only do you gain unique SIM technology. You also maintain the use of standard dyes and fluorescent proteins, the ability to perform simultaneous two-color imaging with clean separation between channels, and the flexibility to choose from a variety of imaging modes to best suit the needs of your samples.




Spheroid stained for mitochondria (MitoTracker Green) and nuclei (NucRed Live 647)

Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Capture entire model organisms and tissue sections


ZEISS Lattice SIM 3 fully leverages SIM Apotome technology, to provide the most outstanding optical sectioning at large fields of view with near-isotropic resolution. ZEISS Lattice SIM 3 is your system of choice for fast imaging of larger volumes, such as 3D model organisms, embryos, organoids, or tissue sections. Whether you work with living or fixed samples, ZEISS Lattice SIM 3 provides access to structured illumination microscopy of multicellular organisms with superior penetration depth.



Spheroid stained for mitochondria (MitoTracker Green) and nuclei (NucRed Live 647)

Acquire super-resolution images as fast and gentle as widefield images


Choose between the standard SIM Apotome imaging mode for the highest available resolution (5 phase images required) or the imaging mode with reduced phases for slightly lower resolution but significantly increased speed and gentleness (only 3 phase images required). Combine SIM Apotome with the Leap mode to significantly speed up super-resolution acquisition. SIM Apotome makes even lossless acquisition possible, meaning for every reconstructed image just one raw image is needed.



Spheroid invading collagen matrix; cells are expressing Lifeact-tdTomato; color-coded depth projection

Go from a large-field overview to the super-resolution details

For large sample experiments, such as imaging whole organisms or simultaneous data collection of a wide area of cells – ZEISS Lattice SIM 3 offers the most advantageous combination of a large field of view and super-resolution imaging. SIM Apotome mode in combination with SIM² image reconstruction enables lateral super-resolution down to 140 nm with superior optical sectioning and sensitivity. Additionally, imaging in Lattice SIM mode with a ZEISS 25x multi-immersion objective and subsequent SIM² processing provides similar lateral resolutions with larger fields of view and more flexible adaptation to the refractive index of your sample.



Cos7 cell expressing ER-mStayGold and stained for mitochondria using MitoTracker Red CMXRos

Your Insight into the Technology Behind It



- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

SIM Apotome: Flexible optical sectioning

Live cell imaging with a widefield system often suffers from out-of-focus blur or background signal. These effects can decrease contrast and resolution. ZEISS Lattice SIM 3 fully leverages the benefits of the SIM Apotome technology, enabling structured illumination microscopy for low-magnification objectives to give you fast and gentle optical sectioning for your multicellular samples.

A grid pattern is used to illuminate and rapidly modulate the fluorescence signals in the focal plane. After acquiring three or five images with different grid positions (phases), these frames are combined into a resulting image which contains only information from the focal plane – your optical section.

The SIM Apotome acquisition mode in combination with the SIM² reconstruction algorithm allows you to further adjust the gentleness of fast live-cell imaging with high contrast and resolution. Or use your new optical sectioning speed to increase productivity when acquiring large sample areas or large volumes at different magnifications.



SIM² Apotome: Comparison of widefield (left) and SIM² Apotome (right) single plane images of Cos-7 cells stained for actin (Phalloidin Alexa Fluor 488, green). Objective: LD LCI Plan-Apochromat 25x/0.8 Imm Corr

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Balance your need for speed and resolution

Higher imaging speeds and decreased light exposures are a constant demand in imaging experiments. At the same time, these acquisition settings affect the resolution of the resulting images, and these parameters must be balanced with the desired outcome. To increase speed and decrease light exposure with SIM techniques, the number of phase images acquired for the reconstruction of one final frame/volume are reduced.

The robustness and flexibility of ZEISS Lattice SIM 3 structured illumination patterns plus the image reconstruction software allow a significant reduction to the number of phase images required for SIM Apotome acquisition mode, and, importantly, this only causes a slight decrease in the resolution of the final images. SIM Apotome acquisition can be operated at 3 phase images per frame, increasing the imaging speed by 66%. The increased speed is also advantageous for fast screening of large sample areas, such as tissue sections.


In combination with the Leap mode, the reduced phase acquisition of SIM Apotome further decreases the number of phase images per final frame. SIM Apotome with reduced phases combined with Leap mode is entirely lossless and provides the same number of phase images and processed frames, enabling the gentlest super-resolution imaging possible.



Live yeast cells expressing a vacuole marker tagged with super-folder GFP. Color-coded depth projection. Cells were imaged for 12 hrs; budding events can be observed. Image courtesy of Chris McDonald, University of York, UK.



SIM Apotome volume tile scan image of Arabidopsis root labeled for the Golgi apparatus; time series was recorded for 35 min; color-coded depth projection. Image courtesy of Peter O'Toole, University of York, UK.



Drosophila embryo stained for Fasciclin II (color-coded depth projection) and HRP (cyan) labelling the nervous system. Image courtesy of Ines Hahn, University of York, UK.

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Lattice SIM:

Your 3D super-resolution technique



ZEISS Lattice SIM 3 also includes the Lattice SIM imaging mode optimized for use with a special 25x multi-immersion objective. The sample area is illuminated with a lattice spot pattern instead of grid lines. The lattice pattern provides higher contrast to allow deeper sample penetration and, in combination with SIM², robust image reconstruction with super-resolution down to 140 nm.

SIM² reconstruction:

Double your SIM resolution

Dual-iterative SIM, or SIM², is a groundbreaking image reconstruction algorithm that increases the resolution and sectioning quality of structured illumination microscopy data. SIM² is compatible with all SIM imaging modes and fully integrated in ZEISS ZEN software.

Unlike conventional reconstruction algorithms, SIM² is a two-step image reconstruction algorithm. First, order combination, denoising, and frequency suppression filtering are performed. All the effects resulting from these digital image manipulations are translated into a digital SIM point spread function (PSF). The subsequent iterative deconvolution uses this PSF. Similar to the advantages of using experimental PSF for deconvolution of hardware-based microscopy data, the SIM² algorithm is superior to conventional one-step image reconstruction methods in terms of resolution, sectioning, and robustness.



Lattice SIM: Comparison of widefield and Lattice SIM images of Cos-7 cells stained for actin (Phalloidin Alexa Fluor 488, magenta), microtubules (anti-beta-tubulin Alexa Fluor 568, yellow) and Paxillin (anti-Paxillin Alexa Fluor 647, cyan). Maximum intensity projections. Objective: 25x /0.8 Imm Corr

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Boost the speed of SIM imaging even further

You can further increase the temporal resolution and productivity for 2D and 3D imaging by using the speed enhancement modes. The Burst mode and the Leap mode are compatible with SIM Apotome as well as Lattice SIM acquisition. Combined with SIM² image reconstruction, they enable you to capture highly dynamic processes at exceptional resolution in all three dimensions. For ZEISS Lattice SIM 3, the combination of SIM Apotome mode with reduced phases and Leap mode allows for super-resolution imaging at widefield speed, i.e., after SIM processing, you get one final super-resolution image per acquired raw image.

2D Burst mode:

Get full temporal information


Burst mode processing uses the rolling window approach to let you observe processes in your living samples at up to 255 fps. Since Burst mode is a post-acquisition step, you have the flexibility to use it with previously acquired data sets. You decide how much temporal resolution is required for your data analysis.




3D Leap mode:

Digital sectioning at a new level

For demanding fast imaging in 3D, the Leap mode acquisition enables you to reduce your imaging time and lower the light exposure on your sample. This works by imaging only every third plane, for three-times higher volume imaging speed and three-times fewer light exposures. ZEN reconstructs the entire volume using a pixel reassignment approach.



SIM Apotome / Burst mode



▶ [Click here to view this video](#)

U2OS cell expressing Golgi derived vesicles (tdTomato, magenta) and Rab5a (mEmerald, green). Objective: 40x / 1.4 Oil

SIM Apotome / Leap mode



▶ [Click here to view this video](#)

U2OS cell expressing EB3-tdTomato, recorded with reduced phases. Objective: 40x / 1.4 Oil

Expand Your Possibilities


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

ZEISS ZEN: A journey through different scales

Biological samples often contain different types of information at different length scales. Collecting low to high resolution data in the same sample not only makes you more productive, but also allows you to interconnect your findings and create more accurate biological models based on your experimental findings. With AI Sample Finder, you automatically detect your whole sample even before starting your experiment – ensuring that you won't miss any relevant areas. The ZEN Connect toolkit enables you to combine different experiments recorded with various acquisition modes or systems – placing your experiments into the spatial context of the whole sample.

ZEISS arivis Pro: Advanced image processing and 3D reconstruction

Use the efficient ZEISS arivis Pro software for visualization and quantification of large 3D and 4D data sets. ZEISS arivis Pro not only renders volume images of almost unlimited size, but also provides advanced image processing tools such as volume fusion, channel shift, conventional and machine learning based segmentation, 3D tracking, and neuron tracing. Visualize your quantitative results within ZEISS arivis Pro or export all data for further analysis. The modular structure of ZEISS arivis Pro flexibly adjusts to your needs for advanced image processing and analysis.



Murine brain expressing the neuronal marker Thy1-eGFP was imaged in SIM Apotome and Lattice SIM modes over a Z stack range of 170 µm. Objective for overview image (left): Plan-Neofluar 10x. Objective for inset (right): LD LCI Plan-Apochromat 25x/0.8 Imm Corr.

The ZEN Connect project combines data sets recorded with 10x SIM Apotome, 25x SIM Apotome, and 25x Lattice SIM. The volume rendering on the right-hand side shows a subset of the 25x Lattice SIM data set. Sample courtesy of Herms Lab (MCN, University of Munich, Germany).


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Super-resolution imaging in Neuroscience

Understanding how neurons respond to damage, disease and metabolic change is critical to how we treat neuronal injury and neurodegenerative disease. The synaptic structure and in particular the active zones where synaptic vesicles are released are key players in signal transmission and proper function of neurons. The imaging of active zones require resolution beyond what can be achieved by standard confocal microscopy.


The lab of Prof. Sean Sweeney investigates a novel mutant that is a regulator of neuronal survival and metabolic responses. Nervous system and synapses are co-labelled with synaptotagmins to observe the general structure of the synapse and distribution of the presynaptic vesicles. Super-resolution microscopy helps identify and quantify the differences in synapse structure and composition of the active zones.



Lower half of *Drosophila* slice stained for nervous system and synapses (Anti-HRP, orange). Objective: Plan-Neofluar 10x/0.3 Air. Image courtesy of Prof. Sean Sweeney, University of York, UK.



Also stained in this section are synaptotagmins (Anti-synaptotagmin, cyan). Objective: LD LCI Plan-Apochromat 25x / 0.8 Imm Corr. The same region of interest was imaged with SIM Apotome and Lattice SIM for comparison. Image courtesy of Prof. Sean Sweeney, University of York, UK.




ZEISS Lattice SIM 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Super-resolution imaging in Immunology

Immunofluorescence of tissue sections is commonly used in immunological research to investigate distribution of and interactions between pathogens and immune cells, all with the aim to develop novel therapies for pathogenic diseases. For compelling results, it is not only crucial to image complete sections as to not miss relevant areas but also to image with enough resolution to identify and quantify individual events.


In the application example shown here, skin tissue sections were imaged with Lattice SIM 3 to investigate distribution of CD8 cells relative to Leishmania parasite infection sites. The enlarged area is a digital zoom-in only, meaning that it is possible to zoom into any region of the overview image and quantify cell nuclei, CD8 cells and Leishmania parasites.




Whole section of skin tissue stained for cell nuclei (cyan) and CD8 cells (yellow). Objective: LD LCI Plan-Apochromat 25x/0.8 Imm Corr. Image courtesy: Helen Ashwin, Department of Biology, University of York, UK.



Region of interests of a skin tissue section stained for cell nuclei (cyan), CD8 cells (yellow) and Leishmania parasites (magenta). Objective: LD LCI Plan-Apochromat 25x/0.8 Imm Corr.



Digital zoom into the image on the left. Parasites can be visualized and quantified in each cell of the section. Image courtesy: Helen Ashwin, Department of Biology, University of York, UK.




ZEISS Lattice SIM 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Super-resolution imaging of living yeast

Live yeast cells are among the most challenging samples in fluorescence microscopy. They are extremely light sensitive and with 4–5 μm in diameter smaller than most used cell lines, e.g., human or mouse. Furthermore, yeast cells grow in suspension; they can move freely in the culture dish and are of spherical shape, without clearly defined orientation. Tackling all these challenges requires extremely gentle and fast imaging combined with high resolution in all spatial dimensions.

SIM² Apotome is the perfect tool to image live yeast cells with super-resolution, yet fast and gentle enough to observe the cells over extended periods of time. The example on the right clearly demonstrates this unique capabilities of SIM² Apotome. Various subcellular compartments (surface marker, endosomes, vacuole, endoplasmic reticulum) were tagged with superfolder GFP and imaged for 12 hours. Yeast cells reproduce quickly, about once every 90 min, by a process called 'budding'. In each of the video sequences, multiple cycles of budding as well as subcellular details and dynamics can be observed.



Multi-well 12-hour time lapse microscopy of live yeast cells expressing superfolder-GFP tagged proteins, color-coded depth projections. Top left: surface protein marker, top right: endosomes, bottom left: vacuole, bottom right: endoplasmic reticulum. Objective: Plan-Apochromat 40 \times /1.4 Oil. Image courtesy of Chris McDonald, University of York, UK.


ZEISS Lattice SIM 3 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Image large volumes with finest detail even at depth

SIM Apotome combined with reduced phases and leap mode allows you to image large volumes extremely fast and efficiently. Crunch through large volumes fast by recording only one raw image per final reconstructed image. Select regions of interest, switch objectives and use Lattice SIM to gain super-resolution images with a lateral resolution down to 140 nm within the context of the whole sample.

A novel clearing and embedding technology developed by Prof. Tang and his team (Hsiao et al., Nature Communications 2023) combined with the advantages of SIM Apotome acquisition and excellent image reconstruction technology enabled us to image an entire mouse intestine section of 3 mm × 4 mm and ~200 µm thickness within a couple of minutes. Networks of blood vessels and nerves can be visualized with finest details even at depth.



▶ [Click here to view this video](#)



▶ [Click here to view this video](#)

Mouse small intestine in A-ha Polymer labeled for blood vessels (Alexa Fluor 488) and nerves (Alexa Fluor 647); anti-fade labeling. Objective: Plan-Neofluar 10×/0.3 Air (left) and LD LCI Plan-Apochromat 25× Imm Corr (right). Sample courtesy of Prof. Shiue-Cheng (Tony) Tang, Institute of Biotechnology & Department of Medical Science, National Tsing Hua University, Taiwan

The Lattice SIM product family

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service

Address all your super-resolution needs across scales

The ZEISS Lattice SIM product family gives you full access to super-resolution imaging for all research areas, from fast optical sectioning to the detection of highly dynamic processes and quantification at the molecular level.



ZEISS Lattice SIM 3

Reveal cellular behavior and inter-cellular dynamics

Lattice SIM 3 is specifically designed to meet the imaging requirements of multicellular organisms and tissue sections. This system exploits the full potential of the SIM Apotome technology: fast optical sectioning at superior quality, large fields of view with access to smaller regions of interest, near-isotropic resolution, and the gentlest super-resolution imaging possible.



ZEISS Lattice SIM 5

Reveal the vibrant sub-organelle network of life

ZEISS Lattice SIM 5 has been optimized for single cell imaging as well as capturing subcellular structures and their dynamics. Powered by the Lattice SIM technology and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 provides you with outstanding super-resolution capabilities down to 60 nm in both living and fixed cells.



ZEISS Elyra 7 with Lattice SIM

Reveal life across scales – down to molecular details

ZEISS Elyra 7 includes several microscopy techniques: Lattice SIM², SIM² Apotome, SMLM and TIRF. You can combine these techniques to multiply the insights from your specimen and to correlate the acquired data. With its focus on single molecule localization microscopy, ZEISS Elyra 7 gives you resolution excellence down to the molecular level.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



1 Microscope

- ZEISS Axio Observer 7 (inverse stand)
- Stage top incubation
- Motorized XY stepper scanning stage
- Z-Piezo stage insert
- 1 camera port for camera or Duolink

2 Objectives

- Plan-Apochromat 40×/1.4 Oil (DIC*)
- C-Apochromat 40×/1.2 W
- LD LCI Plan-Apochromat 25×/0.8 Imm Corr
- Plan-Apochromat 20×/0.8 Air
- EC Plan-Neofluar 10×/0.3 Air

3 Lattice SIM 3 Illumination and Detection

- Fiber coupled diode pumped solid state lasers
- Available lines:
 - 405 nm diode (50 mW),
 - 488 nm diode (50 mW),
 - 561 nm diode (SHG) (50 mW),
 - 640 nm diode (50 mW)
- ZEISS AxioCam 820 CMOS camera


4 Software

- ZEN (blue edition)
- SIM toolkit

* DIC indicates type of objective not imaging modality


System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Incubation Setup

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Microscope

Stand	ZEISS Axio Observer 7, SR for Lattice SIM, motorized inverted microscope for super-resolution microscopy
Z-drive	DC servo motor, opto-electronically coded; smallest Z step: 25 nm
XY stepper scanning stage	Motorized, stepper motor with 2 mm spindle pitch; travel range: 130 mm × 100 mm; max. speed: 50 mm/s; Resolution: 0.1 µm; reproducibility: ± 1 µm; absolute accuracy: ± 5 µm; Suitable for mounting frames K 160 × 110 mm and Z-Piezo stage insert; compatible with objectives' autocorr
Z-Piezo stage insert	For XY scanning stage; max travel range: 500 µm; smallest Z step size: 5 nm; Level-adjustable stage insert for frame inserts (sample holders) and multi-well plates; Sample holders available for 3" × 1" standard slides, LabTek chambers; 35 – 40 mm glass-bottom dishes; Universal stage insert for various carrier formats

Optical Filters

Filter sets reflector turret	Flexible filter set available for simultaneous multi-channel acquisition; Filter set with four precisely mounted ACR-coded filter modules for super-resolution microscopy on a motorized six-position turret; Two positions in the turret compatible with standard Push & Click filter modules, e.g., for visual sample observation
Dual filter set for Duolink	Filter sets are optimized for one color (SOLO), dual color (DUO) and four color (QUAD) applications
Filter slider	Manual filter slider with Bertrand lens; fits into the slit below the objective turret

Lasers

Laser module	Laser coupling with polarization-maintaining single mode fiber (no adjustment of laser coupling by users required)
Laser lines	405 nm (50 mW), 488 nm (50 mW), 561 nm (50 mW), 640 nm (50 mW); 405, 488 & 642 nm: diode lasers (DL); 561 nm: frequency doubled diode laser (FDDL); Direct modulation @ 500:1

Cameras

CMOS	ZEISS AxioCam 820 mono; sensor pixel count: 4512 × 4512 = 20 megapixel, effective: 3072 × 3072; pixel size: 2.74 µm × 2.74 µm; QE: up to 86 % (@460 nm); binning: 1 × 1, 2 × 2 (default), 4 × 4; gain: 1 × (min), 2 ×, 4 × (opt), 8 ×, 16 × (max); active cooling, regulated sensor temperature: 25 °C; bit depth: 14 Bit; frame rate: 28 fps, 75 fps (2 × 2 binning) @ full frame
------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Lattice SIM 3

Illumination module	Illumination module attached to rear port of microscope stand; fully motorized SIM imaging; Two different grating frequencies for SIM Apotome for best match to objective and wavelength; one grating for Lattice SIM; Motorized exchange of gratings in multi-color SIM Apotome; fast piezo actuated phase stepping of gratings.
Camera	Up to two CMOS cameras (ZEISS Axiacam 820) mounted on right side port
Imaging modes	Widefield mode for illumination with thermal light source or LED; laser widefield mode for illumination with laser; Lattice SIM mode using two-dimensional lattice grid; SIM Apotome mode using one-dimensional line grid
Objective lenses (Lattice SIM)	LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*, ACR ⁽¹⁾ coding
Objective lenses (SIM Apotome)	Plan-Apochromat 40×/1.4 Oil; C-Apochromat 40×/1.2 W; LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; Plan-Apochromat 20×/0.8 Air; EC Plan-Neofluar 10×/0.3 Air
Resolution (Lattice SIM/Lattice SIM ²)	Lateral resolution (XY): down to 210 nm/140 nm (typical experimental FWHM values with objective lens LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; with sub-resolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Resolution (SIM/SIM ² Apotome)	Lateral resolution (XY): down to 320/265 nm for 25× (typical experimental FWHM values with sub-resolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Multi-color (Lattice SIM and SIM Apotome)	Detection of up to four different fluorescent labels (sequential detection) and simultaneous dual-color detection with Duolink
Max. field of view (Lattice SIM)	204.3 × 204.3 μm ² , full-frame recording (1536 × 1536 effective px) with LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*
Max. field of view (SIM Apotome)	163.44 × 163.44 μm ² , full frame recording (1536 × 1536 effective px) with Plan-Apochromat 40×/1.4 Oil; 261.51 × 261.51 μm ² , full frame recording with LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; 255.58 × 255.58 μm ² , full frame recording with Plan-Apochromat 20×/0.8 Air; 653.8 × 653.8 μm ² , full frame recording with EC Plan-Neofluar 10×/0.3 Air
Acquisition speed (Lattice SIM)	19 SIM image frames per second at 512 × 512 px resolution and 1 ms exposure time (13 phase images per one SIM image) 28 SIM image frames per second at 512 × 512 px resolution and 1 ms exposure time (9 phase images per one SIM image)
Acquisition speed (SIM Apotome)	51 sectioned frames per second at 512 × 512 px resolution and 1 ms exposure time (camera limited) (5 phase images per one sectioned image); 85 sectioned frames per second at 512 × 512 px resolution and 1 ms exposure time (camera limited) (3 phase images per one sectioned image)
Leap mode and Burst mode	Leap and Burst modes are combinable with both Lattice SIM and SIM Apotome. Leap mode increases the frame rate by a factor of 3 for 3D image acquisition. Max. 255 image frames per second at 512 × 512 px resolution and 1 ms exposure time are available for 2D data after Burst processing.
Data recording and analysis (Lattice SIM and SIM Apotome)	Full software control of SIM imaging; Multi-tracking: sequential multi-channel data acquisition with freely configurable change of gratings (SIM Apotome), or one common grating (Lattice SIM), filters and excitation lasers between tracks; Simultaneous dual-color imaging with one grating; Lattice SIM and SIM Apotome mode imaging in user-defined sub-array regions (ROI imaging); Leap mode for 3 times faster imaging with excellent sectioning; Extension of imaged area possible with tile scanning and stitching; Burst mode processing for 2D time series data sets for Lattice SIM and Apotome mode to increase effective frame rates by a factor of 15 and 5, respectively.

* DIC indicates type of objective not imaging modality

⁽¹⁾ ACR (Automatic Component Recognition); Lattice SIM systems and ZEN imaging software automatically recognize ACR-coded components.

Technical Specifications

-
- › In Brief
-
- › The Advantages
-
- › The Applications
-
- › The System
-
- › **Technology and Details**
-
- › Service
-

Software

Standard	ZEN imaging software (64-bit); operating system: Microsoft Windows 10
	Full software control of image data recording in all imaging modes (including widefield, super-resolution); Software-controlled switching between imaging modes; Full software control of data recording (multi-channel imaging, time series, z-stack); Saving and restoring of user-specific configurations for data recording
SW packages	Required: ZEN module Lattice SIM; ZEN toolkit Advanced Acquisition; ZEN toolkit 3D Optional: ZEN toolkit Deconvolution; ZEN toolkit 2D ; ZEN toolkit Connect; ZEN toolkit AI; ZEN toolkit Developer; Vision package

Accessories

Definite Focus	Holding focus to compensate axial drift, typical z-position accuracy: 30 nm; Specified limits of Definite Focus 3: $0.2 \times \text{DOF}$ (Depth of field: $\text{DOF} \approx \lambda / \text{NA}^2$).
Incubation	Stage top incubation with safety lock
Duolink for attachment of two cameras of the same type	Allows attachment of two cameras of the same type to the microscope.
Storage PC with 81 TByte storage capacity	Direct streaming of data and parallel processing while streaming of data possible



Lattice SIM 3 meets the requirements according to IEC 60825-1:2014 and it is a laser class 1 device.
Interlocks on customer interfaces prevent access to the laser radiation.

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › **Service**

ZEISS Service – Your Partner at All Times

Your microscope system from ZEISS is one of your most important tools. For over 175 years, the ZEISS brand and our experience have stood for reliable equipment with a long life in the field of microscopy. You can count on superior service and support - before and after installation. Our skilled ZEISS service team makes sure that your microscope is always ready for use.

Procurement

- Lab Planning & Construction Site Management
- Site Inspection & Environmental Analysis
- GMP-Qualification IQ/OQ
- Installation & Handover
- IT Integration Support
- Startup Training

Operation


- Predictive Service Remote Monitoring
- Inspection & Preventive Maintenance
 - Software Maintenance Agreements
 - Operation & Application Training
 - Expert Phone & Remote Support
- Protect Service Agreements
 - Metrological Calibration
 - Instrument Relocation
 - Consumables
 - Repairs

New Investment

- Decommissioning
- Trade In


Retrofit

- Customized Engineering
- Upgrades & Modernization
- Customized Workflows via ZEISS arivis Cloud



Please note: Availability of services depends on product line and location

Revealing the vibrant sub-organelle network of life



ZEISS Lattice SIM 5

Your Live Imaging System for Uniform Super-Resolution in All Spatial Dimensions



Seeing beyond

Your Live Imaging System for Uniform Super-Resolution in All Spatial Dimensions

- › In Brief

- › The Advantages

- › The Applications

- › The System

- › Technology and Details

- › Service


The ZEISS Lattice SIM family

Using microscopy to visualize biological structures provides insights into function. When imaging fixed structures, acquisition settings can be optimized for spatial resolution. However, when capturing dynamic events in living samples, higher acquisition speeds and low-light conditions must be balanced with resolution. The ZEISS Lattice SIM family balances sample size, imaging speed, and super-resolution capabilities based on your application – from outstanding optical sectioning of tissues and developing organisms to high-speed imaging of living cells to resolution excellence at the molecular level.

ZEISS Lattice SIM 5

ZEISS Lattice SIM 5 has been optimized for single cell imaging as well as capturing subcellular structures and their dynamics. Powered by the Lattice SIM technology and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 provides you with outstanding super-resolution capabilities down to 60 nm in both living and fixed cells. Additionally, you can choose SIM Apotome imaging mode and a low-magnification objective to achieve fast overview images of your sample before zooming into super-resolution details.

With ZEISS Lattice SIM 5, not only do you gain unique SIM technology. You also maintain the use of standard dyes and fluorescent proteins, the ability to perform simultaneous two-color imaging with clean separation between channels, and the flexibility to choose from a variety of imaging modes to best suit the needs of your samples.




Actin dynamics in a U2OS cell expressing LifeAct-GFP were imaged with the Lattice SIM 3D Leap mode and reduced phases. The movie shows a color-coded depth projection of the volume stack. Objective: Plan-Apochromat 63×/1.4 Oil

Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Capture highly dynamic processes


Equipped with the ZEISS Lattice SIM illumination pattern and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 raises structured illumination microscopy (SIM) to a new level. You will always achieve the best possible results, even when using lower light exposures to protect living specimens. Double the conventional SIM resolution and discriminate the finest subcellular structures that are no more than 60 nm apart. The light-efficient Lattice SIM technology provides the gentlest imaging of living and fixed specimens, giving you not only double spatial resolution compared to classic SIM, but also high temporal resolution with up to 255 fps.



The Lattice SIM² image of Cos-7 cells labeled with phalloidin Alexa Fluor 488 shows the fine structure of the Actin network. Maximum intensity projection of Z stack is shown.

Optimize to the needs of living samples


The flexibility of ZEISS Lattice SIM 5 allows you to balance the needs of your experiment by prioritizing resolution, speed, or by finding the right balance in between. Use the photon budget to enhance lateral resolution well below 100 nm or reduce the number of required raw images to boost acquisition speed and gentleness. ZEISS Lattice SIM 5 has a number of options for reducing raw images which allows you to select for the best acquisition settings that target your desired spatial and temporal resolution.



Lattice SIM²: Time lapse imaging of the endoplasmic reticulum (ER-StayGold) in a Cos-7 cell reveals highly dynamic structural changes. Sample courtesy of Miyawaki Lab, RIKEN Institute, Japan.

Get more reliable experiment results

ZEISS Lattice SIM 5 with SIM² comes with outstanding out-of-focus light suppression, giving you the sharpest sectioning in widefield microscopy, even for highly scattering samples. The SIM² image reconstruction uses a special SIM point spread function to robustly reconstruct all structured-illumination-based acquisition data of your ZEISS Lattice SIM 5 with minimal image artifacts – for both living and fixed samples. Rest assured knowing that you are basing your experimental conclusions on reproducible data generated from a powerful and proven algorithm.



Cos-7 cells stained for microtubules (anti-tubulin Alexa Fluor 488, cyan) and actin (Phalloidin Alexa Fluor 561, orange)

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Lattice SIM:

Your 3D super-resolution technique

In classic SIM, the sample is illuminated with a grid structure that interferes with structures in the sample, creating Moiré fringes. These fringes contain high frequency information – that is, high resolution information – transformed down to low frequencies that can be resolved by the microscope. To achieve this effect in all directions, the sample is imaged at different rotational and translational positions (phases) of the grid pattern. The phase images are deconvolved into the resulting image, which will have twice the resolution in all three dimensions.



In Lattice SIM, the sample is illuminated with a lattice spot pattern instead of grid lines. Due to its intrinsic two-dimensionality, the lattice pattern requires only translational repositioning but no rotation. This leads to a dramatic increase in imaging speed. In addition, the lattice pattern provides higher contrast to allow a more robust image reconstruction. Since the sampling efficiency is doubled compared to classic SIM, half as much light exposure is needed making Lattice SIM a preferred live cell imaging technique.

Lattice SIM

The lattice pattern gives better contrast: you maintain image quality at higher frame rates.



Watch the movie for a quick comparison of classic SIM and Lattice SIM



Lattice SIM: Comparison of widefield and Lattice SIM images of Cos-7 cells stained for actin (Phalloidin Alexa Fluor 568, magenta), microtubules (anti-tubulin Alexa Fluor 488, yellow) and nucleus (Hoechst, blue). Images are maximum intensity projections.
Objective: Plan-Apochromat 63×/1.4 Oil

Your Insight into the Technology Behind It


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

SIM² reconstruction:



Double your SIM resolution

Dual-iterative SIM, or SIM², is a groundbreaking image reconstruction algorithm that increases the resolution and sectioning quality of structured illumination microscopy data. SIM² is compatible with all SIM imaging modes and fully integrated in ZEISS ZEN software.

Unlike conventional reconstruction algorithms, SIM² is a two-step image reconstruction algorithm. First, order combination, denoising, and frequency suppression filtering are performed. All the effects resulting from these digital image manipulations are translated into a digital SIM point spread function (PSF). The subsequent iterative deconvolution uses this PSF. Similar to the advantages of using experimental PSF for deconvolution of hardware-based microscopy data, the SIM² algorithm is superior to conventional one-step image reconstruction methods in terms of resolution, sectioning, and robustness.



GATTA-STED Nanoruler 90B, 80B, 70B and 60B (GATTAquant, Germany) were imaged and processed with Lattice SIM² mode with a 63×/1.4 oil objective. Distances of 90 nm, 80 nm, 70 nm and 60 nm are resolved.



Images of Cos-7 cell stained with anti-alpha-Tubulin Alexa Fluor 488 were processed with the conventional SIM algorithms based on generalized Wiener filter and with the novel SIM² reconstruction. The images show an improvement of resolution for SIM² compared to SIM. The superior sectioning capability of SIM² is shown in the movie. Objective: alpha Plan-Apochromat 100×/1.57 Oil, imaged on Elyra 7 with Lattice SIM.


Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

SIM Apotome: Flexible optical sectioning

Live cell imaging with a widefield system often suffers from out-of-focus blur or background signal. These effects can decrease contrast and resolution. The SIM Apotome acquisition mode uses structured illumination to give you fast optical sectioning of larger volumes with crisp contrast and high resolution in all dimensions. A grid pattern is used to illuminate and rapidly modulate the fluorescence signals in the focal plane. After acquiring three or five images with different grid positions (phases), these frames are combined into a resulting image which contains only information from the focal plane – your optical section.



The SIM Apotome acquisition mode in combination with the SIM² reconstruction algorithm allows you to further tune the gentleness of fast live-cell imaging with high contrast and resolution. Or use your new optical sectioning speed to increase productivity when acquiring large sample areas or large volumes at different magnifications.



Widefield image with out-of-focus light. Signal from the focal plane is encircled by a white dashed line.

SIM Apotome acquisition at 3 or 5 different grid positions

Reconstructed optically sectioned image



SIM² Apotome: Comparison of widefield and SIM² Apotome single plane images of U2OS cells stained for actin (phalloidin Alexa Fluor 647, red), microtubules (anti-alpha-tubulin Alexa Fluor 488, green) and nuclei (Hoechst, blue).
Objective: LD LCI Plan-Apochromat 25×/0.8 Imm Corr

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Balance your need for speed and resolution

Higher imaging speeds and decreased light exposures are a constant demand in imaging experiments. At the same time, these acquisition settings affect the resolution of the resulting images, and these parameters must be balanced with the desired outcome. To increase speed and decrease light exposure with SIM techniques, the number of phase images acquired for the reconstruction of one final frame/volume are reduced.


The robustness and flexibility of ZEISS Lattice SIM 5 structured illumination patterns plus the image reconstruction software allow a significant reduction to the number of phase images required for Lattice SIM acquisition mode, and,

importantly, this only causes a slight decrease in the resolution of the final images. Lattice SIM 5 acquisition can be operated at 9 phase images per frame instead of 13, increasing the imaging speed by 44 %. The increased speed is particularly advantageous for gentle imaging of highly dynamic live cells where slower acquisition would result in motion blur and reduced resolution.


In combination with Leap mode, the reduced phase acquisition of Lattice SIM decreases the number of phase images per final frame. This means that for every reconstructed image, just one raw image is needed, enabling gentle super-resolution imaging that is unprecedented.



Cos-7 cell expressing TMM20-mEmerald (cyan) and EB3-tdTomato (orange) shows dynamic movement of mitochondria and microtubules. Imaged with Lattice SIM. Objective: Plan-Apochromat 63x/1.4 Oil



Actin dynamics of U2OS cells expressing LifeAct-tdTomato (orange) imaged with the Lattice SIM² mode with reduced phases. The cells were embedded in a collagen matrix stained with FastGreen dye (cyan). The image shows a maximum intensity projection of the volume stack. Objective: Plan-Apochromat 63x/1.4 Oil.



Actin dynamics in a U2OS cell expressing LifeAct-GFP were imaged with the SIM Apotome 3D Leap mode and reduced phases. The image shows a maximum intensity projection of the volume stack. Objective: Plan-Apochromat 40x/1.4 Oil

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service



Boost the speed of SIM imaging even further

Lattice SIM 5 already provides you with fast imaging speeds. But you can further increase the temporal resolution and productivity for 2D and 3D imaging by using the speed enhancement modes. The Burst mode and the Leap mode are compatible with Lattice SIM as well as SIM Apotome acquisition. Combined with SIM² image reconstruction, they enable you to capture highly dynamic processes at exceptional resolution in all three dimensions.

2D Burst mode:

Get full temporal information

Burst mode processing uses the rolling window approach to let you observe processes in your living samples at up to 255 fps. Since Burst mode is a post-acquisition step, you have the flexibility to use it with previously acquired data sets. You decide how much temporal resolution is required for your data analysis.





U2OS cell expressing Rab5-mEmerald (green) and tdTomato tagged Golgi associated transport marker (magenta). Simultaneous dual-color acquisition with an exposure time of 1.5 ms/phase for a FOV of 1024×1024 pixel (64 μm × 64 μm). Objective: Plan-Apochromat 63×/1.4 Oil

3D Leap mode:

Digital sectioning at a new level

For demanding fast imaging in 3D, the Leap mode acquisition enables you to reduce your imaging time and lower the light exposure on your sample. This works by imaging only every third plane, for three-times higher volume imaging speed and three-times fewer light exposures. ZEN reconstructs the entire volume using a pixel reassignment approach.



U2OS cell expressing calreticulin-tdTomato to visualize the endoplasmic reticulum. The time series shows a maximum intensity projection of the volume data set. Objective: Plan-Apochromat 63×/1.4 Oil

Expand Your Possibilities


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

ZEISS ZEN: A journey through different scales

Biological samples often contain different types of information at different length scales. Collecting low to high resolution data in the same sample not only makes you more productive, but also allows you to interconnect your findings and create more accurate biological models based on your experimental findings. With AI Sample Finder, you automatically detect your whole sample even before starting your experiment – ensuring that you won't miss any relevant areas. The ZEN Connect toolkit enables you to combine different experiments recorded with various acquisition modes or systems – placing your experiments into the spatial context of the whole sample.

ZEISS arivis Pro: Advanced image processing and 3D reconstruction

Use the efficient ZEISS arivis Pro software for visualization and quantification of large 3D and 4D data sets. ZEISS arivis Pro not only renders volume images of almost unlimited size, but also provides advanced image processing tools such as volume fusion, channel shift, conventional and machine learning based segmentation, 3D tracking, and neuron tracing. Visualize your quantitative results within ZEISS arivis Pro or export all data for further analysis. The modular structure of ZEISS arivis Pro flexibly adjusts to your needs for advanced image processing and analysis.



Murine brain expressing the neuronal marker Thy1-eGFP was imaged in SIM Apotome and Lattice SIM modes over a Z stack range of 170 μm . Objective for overview image (left): Plan-Neofluar 10x. The ZEN Connect project combines data sets recorded with 10x SIM Apotome, 25x SIM Apotome, 40x SIM Apotome and 63x Lattice SIM. The volume rendering on the right-hand side shows a subset of the 63x Lattice SIM data set. Objective: Plan-Apochromat 63x/1.4 Oil. Sample courtesy of Herms Lab (MCN, University of Munich, Germany).

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Simultaneous two-color imaging

Investigation of living samples very often focuses on interactions of different proteins or organelles. Simultaneous imaging of the involved structures is key to proper understanding of these highly dynamic processes. ZEISS Lattice SIM 5 can be equipped with a Duolink adapter to operate two sCMOS cameras in parallel and perform true simultaneous two-color imaging within your entire field of view.

ZEISS Axiocam 820 mono


For excellent performance at a cost-efficient price point, select the ZEISS Axiocam 820 mono camera, which features back-illuminated CMOS sensor with a peak quantum efficiency of 86 %. In combination with its low readout noise, this camera is the perfect choice for imaging faint fluorescence signals in living or fixed samples. The USB 3.0 interface enables high acquisition speeds and exposure times down to 1 ms.



ZEISS Lattice SIM 5 equipped with two ZEISS Axiocam 820 mono cameras

Hamamatsu ORCA-Fusion BT

For outstanding performance, choose the Hamamatsu ORCA-Fusion BT camera. This camera features a scientific CMOS (sCMOS) with a back-thinned sensor enabling peak quantum efficiency of ~95 %. With its ultra-low, uniform readout noise and CoaXPress interface for high acquisition speeds and exposure times down to 1 ms, it yields unparalleled digital imaging results.



ZEISS Lattice SIM 5 equipped with two Hamamatsu ORCA-Fusion BT cameras


ZEISS Lattice SIM 5 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Straightforward multi-color super-resolution imaging

Studying multiprotein complexes requires super-resolution imaging with multiple colors, which is not often attainable with conventional techniques. Lattice SIM² enables you to perform multi-color imaging at resolution down to 60 nm for conventionally stained samples.

The synaptonemal complex is a well-known structure in the nucleus of meiotic cells and consists of two lateral elements, which are connected to a central element by transverse filaments. Due to its small size, three-color imaging of the synaptonemal complex has previously been possible only using complex methods with elaborate sample preparation, such as super-resolution imaging of three-fold expanded samples using the expansion microscopy technique. Lattice SIM² resolves the two strands of SYCP3 (lateral elements) as well as SYCP1-C (C-terminus of transverse filaments) without special sample treatment or staining for distances well below 100 nm. More importantly, the three-color image provides structural information for the distances between the proteins SYCP3 and SYCP1. Even within the SYCP1 protein, the differently labeled N- and C-Terminus can be clearly separated with less than 50 nm resolution between the two labels.



Architecture of threefold labeled synaptonemal complexes from mouse testis visualized via immunolabeling of SYCP3 with SeTau647, SYCP1-C with Alexa Fluor 488 and SYCP1-N with Alexa Fluor 568 and Lattice SIM² mode. Objective: Plan-Apochromat 63x/1.4 Oil. Sample courtesy: Marie-Christin Spindler, AG Prof Ricardo Benavente, Biocenter of the University of Würzburg.

ZEISS Lattice SIM 5 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Observe life's finest details

Understanding biological processes requires observation of living cells or organisms at low light dosage and high spatiotemporal resolution. ZEISS Lattice SIM 5 is your super-resolution system designed for imaging live specimens. Due to its unique lattice structural illumination, it combines high speed imaging with incredible light efficiency, low photon dosage and sensitivity. You can observe cellular, subcellular, and even


sub-organelle structures in living specimens in 2D and 3D over time. Whether you are interested in cytoskeletal dynamics, mitochondrial fusion and fission or budding of the endoplasmic reticulum, Lattice SIM 5 provides you with the necessary live cell compatibility at super-resolution.

Mitochondria are powerhouses of our cells, generating energy in the form of ATP to sustain

the cell. They are highly dynamic organelles that constantly undergo fusion and fission events to ensure proper distribution of ATP across the cell. In order to do their job, they are known to interact with many other subcellular compartments including microtubules, which they ride on to get to their destinations, or the ER, which wraps around mitochondria to initially constrict their diameter before fission events.



U2OS cell expressing Tomm20-mEmerald. Image shows a color-coded projection of the Lattice SIM² volume data set. Objective: Plan-Apochromat 63×/1.4 Oil



Cos-7 cell expressing the endoplasmic reticulum marker Calreticulin-tdTomato (magenta) and mitochondrial marker Tomm20-mEmerald (green) was simultaneously imaged for two colors. The movie shows high dynamic interactions of the ER and mitochondria. Objective: Plan-Apochromat 63×/1.4 Oil


ZEISS Lattice SIM 5 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Observe life's finest details

Study of the components of the cytoskeleton is a prominent research field in biology. Due to the fine structures of these components, for example actin and ER networks or microtubular filaments, imaging far below 100 nm is often performed with super-resolution techniques. Lattice SIM² allows you to gain much more structural information from your samples compared to conventional SIM techniques. It not only operates with a resolution of down to 60 nm

but also provides markedly improves the sectioning quality in your images. This robust image reconstruction method efficiently separates signal and background without the need for tailor-made staining protocols or expert knowledge of complex microscopy techniques. You can take advantage of the easy-to-use Lattice SIM² technology to unveil complex structural information and gain more insights from your experiments.



Actin dynamics in a U2OS cell expressing LifeAct-GFP were imaged with the Lattice SIM 3D Leap mode and reduced phases. The movie shows a color-coded depth projection of the volume stack. Objective: Plan-Apochromat 63×/1.4 Oil



Simultaneous imaging of the endoplasmic reticulum (Calreticulin-tdTomato, magenta) and microtubules (EMTB-3xGFP, green) in a Cos-7 cell reveals highly dynamic interaction of these organelles. Objective: Plan-Apochromat 63×/1.4 Oil

ZEISS Lattice SIM 5 at Work


- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Resolve the details hiding in the depth

The Lattice SIM illumination pattern exhibits both higher contrast and deeper sample penetration as compared to classical SIM. Achieve super-resolution images along with high-quality sectioning even in thick or scattering samples.

A novel clearing and embedding technology developed by Prof. Tang and his team (Hsiao et al., Nature Communications 2023) combined with the robust Lattice SIM illumination pattern and excellent image reconstruction technology

enabled imaging throughout an entire mouse intestine section of ~200 μm thickness. Networks of blood vessels and nerves can be visualized with finest details even at this depth.



Mouse small intestine in A-ha Polymer labeled for blood vessels (Alexa Fluor 488) and nerves (Alexa Fluor 647); anti-fade labeling. Left: Overview image of the whole section recorded with SIM Apotome, blood vessels: color-coded depth projection, nerves: cyan. Objective: Plan-Neofluar 10 \times /0.3 Air. Top right: Digital zoom into overview image. Objective: Plan-Neofluar 10 \times /0.3 Air. Bottom right: Selected region of interest imaged with Lattice SIM, blood vessels: orange. Objective: Plan Apochromat 63 \times /1.4 Oil. Sample courtesy of Prof. Shiue-Cheng (Tony) Tang, Institute of Biotechnology & Department of Medical Science, National Tsing Hua University, Taiwan

The Lattice SIM product family

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service

Address all your super-resolution needs across scales

The ZEISS Lattice SIM product family gives you full access to super-resolution imaging for all research areas, from fast optical sectioning to the detection of highly dynamic processes and quantification at the molecular level.



ZEISS Lattice SIM 3

Reveal cellular behavior and inter-cellular dynamics

Lattice SIM 3 is specifically designed to meet the imaging requirements of multicellular organisms and tissue sections. This system exploits the full potential of the SIM Apotome technology: fast optical sectioning at superior quality, large fields of view with access to smaller regions of interest, near-isotropic resolution, and the gentlest super-resolution imaging possible.



ZEISS Lattice SIM 5

Reveal the vibrant sub-organelle network of life

ZEISS Lattice SIM 5 has been optimized for single cell imaging as well as capturing subcellular structures and their dynamics. Powered by the Lattice SIM technology and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 provides you with outstanding super-resolution capabilities down to 60 nm in both living and fixed cells.




ZEISS Elyra 7 with Lattice SIM

Reveal life across scales – down to molecular details

ZEISS Elyra 7 includes several microscopy techniques: Lattice SIM², SIM² Apotome, SMLM and TIRF. You can combine these techniques to multiply the insights from your specimen and to correlate the acquired data. With its focus on single molecule localization microscopy, ZEISS Elyra 7 gives you resolution excellence down to the molecular level.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



1 Microscope

- ZEISS Axio Observer 7 (inverse stand)
- Stage top incubation
- Motorized XY stepper scanning stage
- Z-Piezo stage insert
- 1 camera port for camera or Duolink

2 Objectives

- alpha Plan-Apochromat 63×/1.46 Oil
- Plan-Apochromat 63×/1.4 Oil (DIC*)
- C-Apochromat 63×/1.2 Water (DIC*)
- Plan-Apochromat 40×/1.4 Oil (DIC*)
- C-Apochromat 40×/1.2 W
- LD LCI Plan-Apochromat 25×/0.8 Imm Corr
- Plan-Apochromat 20×/0.8 Air
- EC Plan-Neofluar 10×/0.3 Air

3 Lattice SIM 5 Illumination and Detection

- Fiber-coupled diode pumped solid state lasers
- Available lines:
 - 405 nm diode (50 mW),
 - 488 nm diode (50 mW),
 - 561 nm diode (SHG) (50 mW),
 - 640 nm diode (50 mW)
- ZEISS AxioCam 820 mono CMOS camera
- Hamamatsu ORCA-Fusion BT sCMOS camera


4 Software

- ZEN (blue edition)
- SIM toolkit

* DIC indicates type of objective not imaging modality


System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Incubation Setup

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Technical Specifications

-
- › In Brief
-
- › The Advantages
-
- › The Applications
-
- › The System
-
- › **Technology and Details**
-
- › Service
-

Microscope

Stand	ZEISS Axio Observer 7, motorized inverted microscope for super-resolution microscopy
Z-drive	DC servo motor, opto-electronically coded; smallest Z step: 25 nm
XY stepper scanning stage	Motorized, stepper motor with 2 mm spindle pitch; travel range: 130 mm × 100 mm; max. speed: 50 mm/s; Resolution: 0.1 µm; reproducibility: ± 1 µm; absolute accuracy: ± 5 µm; Suitable for mounting frames K 160 × 110 mm and Z-Piezo stage insert; compatible with objectives' autocorr
Z-Piezo stage insert	For XY scanning stage, max travel range: 500 µm; smallest Z step size: 5 nm; Level-adjustable stage insert for frame inserts (sample holders) and multi-well plates; Sample holders available for 3" × 1" standard slides, LabTek chambers; 35–40 mm glass-bottom dishes; Universal stage insert for various carrier formats

Optical Filters

Filter sets reflector turret	Flexible filter set available for simultaneous multi-channel acquisition Filter set with four precisely mounted ACR-coded filter modules for super-resolution microscopy on a motorized six-position turret; Two positions in the turret compatible with standard Push & Click filter modules, e.g., for visual sample observation
Dual filter set for Duolink	Filter sets are optimized for one color (SOLO), dual color (DUO) and four color (QUAD) applications
Filter slider	Manual filter slider with Bertrand lens; fits into the slit below the objective turret

Lasers

Laser module for Lattice SIM 5	Laser coupling with polarization-maintaining single mode fiber (no adjustment of laser coupling by users required)
Laser lines	405 nm (50 mW), 488 nm (50 mW), 561 nm (50 MW), 640 nm (50 mW); 405, 488 & 640 nm: diode lasers (DL); 561 nm: frequency doubled diode laser (FDDL); Direct modulation @ 500:1

Cameras

CMOS	ZEISS AxioCam 820 mono; sensor pixel count: 4512 × 4512 = 20 megapixel, effective: 3072 × 3072; pixel size: 2.74 µm × 2.74 µm; QE: up to 86 % (@460 nm); binning: 1 × 1, 2 × 2 (default), 4 × 4; gain: 1 × (min), 2 ×, 4 × (opt), 8 ×, 16 × (max); active cooling, regulated sensor temperature: 25 °C; bit depth: 14 Bit; frame rate: 28 fps, 75 fps (2 × 2 binning) @ full frame
sCMOS	Hamamatsu ORCA-Fusion BT; sensor pixel count: 2304 × 2304, effective: 1304 × 1304; pixel size 6.5 µm × 6.5 µm; QE: up to 95 % @ 540 nm; water cooling (regulated sensor temperature –8 °C); dynamic range: 16 bit; binning: 1 × 1, 2 × 2, 4 × 4; Frame rates: 89 fps (fast) @ full frame Liquid cooling system for sCMOS camera (Hamamatsu ORCA-Fusion BT)

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Lattice SIM 5

Illumination module	Illumination module attached to rear port of microscope stand; fully motorized SIM imaging; Five different grating frequencies for Lattice SIM for optimal matching of illumination pattern to laser wavelength and objective lens; Motorized exchange of gratings in multi-color Lattice SIM; one grating for SIM Apotome; fast piezo actuated phase stepping of gratings
Camera	Up to two CMOS (ZEISS Axiocam 820) or sCMOS (Hamamatsu ORCA-Fusion BT) cameras mounted on right side port
Imaging modes	Widefield mode for illumination with thermal light source or LED; laser widefield mode for illumination with laser; Lattice SIM mode using two-dimensional lattice grid; SIM Apotome mode using one-dimensional line grid
Objective lenses (Lattice SIM)	Plan-Apochromat 63×/1.40 Oil DIC*; C-Apochromat 63×/1.20 W Corr; alpha Plan-Apochromat 63×/1.46 Oil, ACR ⁽¹⁾ coding
Objective lenses (SIM Apotome)	Plan-Apochromat 40×/1.4 Oil; C-Apochromat 40×/1.2 W; LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; Plan-Apochromat 20×/0.8 Air; EC Plan-Neofluar 10×/0.3 Air
Resolution (Lattice SIM/Lattice SIM ²)	Lateral resolution (XY): down to 120/60 nm (typical experimental FWHM values with objective lens Plan-Apochromat 63×/1.40 Oil DIC*, sub-resolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Resolution (SIM/SIM ² Apotome)	Lateral resolution (XY): down to 320/265 nm for 25× (typical experimental FWHM values with sub-resolution beads of 100 nm diameter and excitation at 488 nm)
Multi-color (Lattice SIM and SIM Apotome)	Detection of up to four different fluorescent labels (sequential detection) and simultaneous dual-color detection with Duolink
Max. field of view (Lattice SIM) @ ORCA-Fusion BT	103.21 × 103.21 μm ² , full-frame recording (1288 × 1288 effective px) with Plan-Apochromat 63×/1.40 Oil DIC*
Max. field of view (SIM Apotome) @ ORCA-Fusion BT	127 × 127 μm ² , full frame recording (1288 × 1288 effective px) with Plan-Apochromat 40×/1.40 Oil; 203.20 × 203.20 μm ² , full frame recording with LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; 254 × 254 μm ² , full frame recording with Plan-Apochromat 20×/0.8 Air; 651 × 651 μm ² , full frame recording with EC Plan-Neofluar 10×/0.3 Air
Acquisition speed (Lattice SIM)	19 SIM image frames per second at 512 × 512 px resolution and 1 ms exposure time (13 phase images per one SIM image) 28 SIM image frames per second at 512 × 512 px resolution and 1 ms exposure time (9 phase images per one SIM image)
Acquisition speed (SIM Apotome)	51 sectioned frames per second at 512 × 512 px resolution and 1 ms exposure time (camera limited) (5 phase images per one sectioned image); 85 sectioned frames per second at 512 × 512 px resolution and 1 ms exposure time (camera limited) (3 phase images per one sectioned image)
Leap mode and Burst mode	Leap and Burst modes are combinable with both the Lattice SIM and SIM Apotome. Leap mode increases the frame rate by a factor of 3 for 3D image acquisition. Max. 255 image frames per second at 512 × 512 px resolution and 1 ms exposure time are available for 2D data after Burst processing.
Data recording and analysis (Lattice SIM and SIM Apotome)	Full software control of Lattice SIM imaging; Multi-tracking: sequential multi-channel data acquisition with freely configurable change of gratings (Lattice SIM), or one common grating (SIM Apotome mode), filters and excitation lasers between tracks; Simultaneous dual-color imaging with one grating; Lattice SIM and SIM Apotome mode imaging in user-defined sub-array regions (ROI imaging); Leap mode for 3 times faster imaging with excellent sectioning; Extension of imaged area possible with tile scanning and stitching; Burst mode processing for 2D time series data sets for Lattice SIM and Apotome mode to increase effective frame rates by a factor of 15 and 5, respectively.

* DIC indicates type of objective not imaging modality

⁽¹⁾ ACR (Automatic Component Recognition); Lattice SIM systems and ZEN imaging software automatically recognize ACR-coded components.

Technical Specifications

-
- › In Brief
-
- › The Advantages
-
- › The Applications
-
- › The System
-
- › **Technology and Details**
-
- › Service
-

Software

Standard	ZEN imaging software (64-bit); operating system: Microsoft Windows 10
	Full software control of image data recording in all imaging modes (including widefield, super-resolution); Software-controlled switching between imaging modes; Full software control of data recording (multi-channel imaging, time series, z-stack); Saving and restoring of user-specific configurations for data recording
SW packages	Required: ZEN package; ZEN module Lattice SIM; ZEN toolkit Advanced Acquisition; ZEN toolkit 3D, ZEN toolkit 2D Optional: ZEN toolkit Deconvolution; ZEN toolkit Connect; ZEN toolkit AI; ZEN toolkit Developer; Vision package

Accessories

Definite Focus	Holding focus to compensate axial drift, typical z-position accuracy: 30 nm; Specified limits of Definite Focus 3: $0.2 \times \text{DOF}$ (Depth of field: $\text{DOF} \approx \lambda / \text{NA}^2$).
Incubation	Stage top incubation with safety lock
Duolink for attachment of two cameras of the same type	Allows attachment of two cameras of the same type to the microscope.
Storage PC with 81 TByte storage capacity	Direct streaming of data and parallel processing while streaming of data possible



Lattice SIM 5 meets the requirements according to IEC 60825-1:2014 and it is a laser class 1 device.
Interlocks on customer interfaces prevent access to the laser radiation.

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › Technology and Details
- › **Service**

ZEISS Service – Your Partner at All Times

Your microscope system from ZEISS is one of your most important tools. For over 175 years, the ZEISS brand and our experience have stood for reliable equipment with a long life in the field of microscopy. You can count on superior service and support - before and after installation. Our skilled ZEISS service team makes sure that your microscope is always ready for use.

Procurement

- Lab Planning & Construction Site Management
- Site Inspection & Environmental Analysis
- GMP-Qualification IQ/OQ
- Installation & Handover
- IT Integration Support
- Startup Training

Operation


- Predictive Service Remote Monitoring
- Inspection & Preventive Maintenance
 - Software Maintenance Agreements
 - Operation & Application Training
 - Expert Phone & Remote Support
- Protect Service Agreements
 - Metrological Calibration
 - Instrument Relocation
 - Consumables
 - Repairs

New Investment

- Decommissioning
- Trade In

Retrofit

- Customized Engineering
- Upgrades & Modernization
- Customized Workflows via ZEISS arivis Cloud



Please note: Availability of services depends on product line and location

Revealing life across scales – down to molecular details



ZEISS Elyra 7 with Lattice SIM

Your Complete Super-Resolution System with Unprecedented Resolution



Seeing beyond

Your Complete Super-Resolution System with Unprecedented Resolution

- › In Brief

- › The Advantages

- › The Applications

- › The System

- › Technology and Details

- › Service

The ZEISS Lattice SIM family

Using microscopy to visualize biological structures provides insights into function. When imaging fixed structures, acquisition settings can be optimized for spatial resolution. However, when capturing dynamic events in living samples, higher acquisition speeds and low-light conditions must be balanced with resolution. The ZEISS Lattice SIM family balances sample size, imaging speed, and super-resolution capabilities based on your application – from outstanding optical sectioning of tissues and developing organisms to high-speed imaging of living cells to resolution excellence at the molecular level.

ZEISS Elyra 7 with Lattice SIM

ZEISS Elyra 7 includes a wealth of microscopy techniques to meet your experimental needs across scales, optimally matching resolution, speed, and sensitivity requirements to your demanding samples. Employ SIM Apotome for fast optical sectioning, Lattice SIM for super-resolution imaging, SIM² image reconstruction for resolution excellence down to 60 nm, as well as SMLM and TIRF for investigations at the molecular level. You can combine these techniques to multiply the insights from your specimen and to correlate the acquired data.

With ZEISS Elyra 7, not only do you gain unique SIM technology. You also maintain the use of standard dyes and fluorescent proteins, the ability to perform simultaneous two-color imaging with clean separation between channels, and the flexibility to choose from a variety of imaging modes to best suit the needs of your samples.




Simpler. More Intelligent. More Integrated.

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Resolution excellence with Lattice SIM²

With SIM², a novel image reconstruction algorithm raising the SIM technology to a new level, you make the best use of the available photon budget. You can now double the conventional SIM resolution and discriminate the finest subcellular structures, even those no more than 60 nm apart. Lattice SIM² comes with outstanding out-of-focus light suppression, giving you the sharpest sectioning in wide-field microscopy even for highly scattering samples. SIM² image reconstruction robustly reconstructs all structured-illumination-based acquisition data of your Elyra 7 – with minimal artefacts – for living and fixed samples.




The Lattice SIM² image of Cos-7 cells labeled with phalloidin Alexa 488 shows the fine structure of the Actin network. Maximum intensity projection of Z stack is shown.

Speed and efficiency for your experiments

While doubling the classic SIM resolution, the light-efficient Lattice SIM² gives you gentle imaging of living and fixed specimens at high speeds of up to 255 fps.

Combine Lattice SIM² and SIM² Apotome with Burst and Leap modes to make super-resolution acquisition faster than ever before. With SIM Apotome mode, even lossless acquisition can be achieved, meaning for every reconstructed image just one raw image is needed!

Or make use of Elyra 7 Duolink to image two differently stained structures simultaneously and use the multiple colors to boost resolution even further.




Lattice SIM²: Time lapse imaging of the endoplasmic reticulum (Calreticulin-tdTomato) in a Cos-7 cell reveals highly dynamic structural changes.

Flexibility for your research

Elyra 7 handles virtually all types of samples, including photo-sensitive cell cultures, scattering *C. elegans* and plants or tissue sections of up to 100 µm thickness.

Elyra 7 includes several microscopy techniques: Lattice SIM², SIM² Apotome, widefield, SMLM and TIRF. You can correlate images of the same sample acquired using any or all of all these techniques to multiply the insights from your specimen.

You can even combine Elyra 7 with a variety of other imaging systems such as LSM with Airyscan or scanning electron microscopy in a correlative workflow.



SMLM: *Xenopus laevis* A6 cells (epithelial kidney cells). Gp120, a nuclear pore complex protein arranged with eightfold symmetry was labeled with Alexa Fluor 647.

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Lattice SIM:

Your 3D super-resolution technique

In classic SIM, the sample is illuminated with a grid structure that interferes with structures in the sample, creating Moiré fringes. These fringes contain high frequency information – that is, high resolution information – transformed down to low frequencies that can be resolved by the microscope. To achieve this effect in all directions, the sample is imaged at different rotational and translational positions (phases) of the grid pattern. The phase images are deconvolved into the resulting image, which will have twice the resolution in all three dimensions.



In Lattice SIM, the sample is illuminated with a lattice spot pattern instead of grid lines. Due to its intrinsic two-dimensionality, the lattice pattern requires only translational repositioning but no rotation. This leads to a dramatic increase in imaging speed. In addition, the lattice pattern provides higher contrast to allow a more robust image reconstruction. Since the sampling efficiency is doubled compared to classic SIM, half as much light exposure is needed making Lattice SIM a preferred live cell imaging technique.

Lattice SIM

The lattice pattern gives better contrast: you maintain image quality at higher frame rates.



Watch the movie for a quick comparison of classic SIM and Lattice SIM



Lattice SIM: Comparison of widefield and Lattice SIM images of a Cos-7 cell undergoing mitosis stained for actin (Phalloidin Alexa Fluor 568, magenta), microtubules (anti-beta-tubulin Alexa Fluor 488, yellow) and nucleus (Hoechst, blue). Images are maximum intensity projections of 30 planes of a total depth of 3.19 μm . Objective: Plan-Apochromat 63 \times /1.4 Oil

Your Insight into the Technology Behind It


- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

SIM² reconstruction:


Double your SIM resolution

Dual-iterative SIM, or SIM², is a groundbreaking image reconstruction algorithm that increases the resolution and sectioning quality of structured illumination microscopy data. SIM² is compatible with all SIM imaging modes and fully integrated in ZEISS ZEN software.

Unlike conventional reconstruction algorithms, SIM² is a two-step image reconstruction algorithm. First, order combination, denoising, and frequency suppression filtering are performed. All the effects resulting from these digital image manipulations are translated into a digital SIM point spread function (PSF). The subsequent iterative deconvolution uses this PSF. Similar to the advantages of using experimental PSF for deconvolution of hardware-based microscopy data, the SIM² algorithm is superior to conventional one-step image reconstruction methods in terms of resolution, sectioning, and robustness.



GATTA-STED Nanoruler 90B, 80B, 70B and 60B (GATTAquant, Germany) were imaged and processed with Lattice SIM² mode with a 63×/1.4 oil objective. Distances of 90 nm, 80 nm, 70 nm and 60 nm are resolved.



Images of Cos-7 cell stained with anti-alpha-Tubulin Alexa Fluor 488 were processed with the conventional SIM algorithms based on generalized Wiener filter and with the novel SIM² reconstruction. The images show an improvement of resolution for SIM² compared to SIM. The superior sectioning capability of SIM² is shown in the movie. Objective: Plan-Apochromat 63× / 1.4 Oil


Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


SIM Apotome: Flexible optical sectioning

Live cell imaging with a widefield system often suffers from out-of-focus blur or background signal. These effects can decrease contrast and resolution. The SIM Apotome acquisition mode uses structured illumination to give you fast optical sectioning of larger volumes with crisp contrast and high resolution in all dimensions. A grid pattern is used to illuminate and rapidly modulate the fluorescence signals in the focal plane. After acquiring three or five images with different grid positions (phases), these frames are combined into a resulting image which contains only information from the focal plane – your optical section.


The SIM Apotome acquisition mode in combination with the SIM² reconstruction algorithm allows you to further tune the gentleness of fast live-cell imaging with high contrast and resolution. Or use your new optical sectioning speed to increase productivity when acquiring large sample areas or large volumes at different magnifications.





Widefield image with out-of-focus light. Signal from the focal plane is encircled by a white dashed line.



SIM Apotome acquisition at 3 or 5 different grid positions



Reconstructed optically sectioned image



SIM² Apotome: Comparison of widefield and SIM² Apotome single plane images of Cos-7 cells stained for microtubules (anti- α -tubulin Alexa Fluor 488, green) and nuclei (Hoechst, blue). Objective: LD LCI Plan-Apochromat 25 \times /0.8 Imm Corr


Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Balance your need for speed and resolution

Higher imaging speeds and decreased light exposures are a constant demand in imaging experiments. At the same time, these acquisition settings affect the resolution of the resulting images, and these parameters must be balanced with the desired outcome. To increase speed and decrease light exposure with SIM techniques, the number of phase images acquired for the reconstruction of one final frame/volume are reduced. The robustness and flexibility of ZEISS Elyra 7 structured illumination patterns plus the image reconstruction software allow a significant reduction to the number of phase images required for Lattice SIM acquisition mode, and, importantly, this only causes a slight decrease in the resolution of the final images. Elyra 7 Lattice SIM acquisition


can be operated at 9 phase images per frame instead of 13; SIM Apotome acquisition can be operated at 3 phase images per frame instead of 5; increasing the imaging speed by 44 % or 66 % respectively. The increased speed is particularly advantageous for gentle imaging of highly dynamic live cells where slower acquisition would result in motion blur and reduced resolution. In combination with Leap mode, the reduced phase acquisition of Lattice SIM decreases the number of phase images per final frame. This means that you only need three times as many phase images as the resulting number of reconstructed full frames, enabling gentle super-resolution imaging that is unprecedented.



Cos-7 cell expressing EMTB-3xGFP (green) and EB3-tdTomato (magenta) shows dynamic movement of microtubules. Imaged in Lattice SIM 9 phase mode. Objective: Plan-Apochromat 63x/1.4 Oil



Actin dynamics in a Cos-7 cell expressing LifeAct-tdTomato were imaged with the SIM Apotome 3D Leap mode over time. The image shows a maximum intensity projection of 30 planes over 3.4 μm depth. Objective: Plan-Apochromat 40x/1.4 Oil



SIM Apotome volume tile scan image of a Helianthus section. Sample: "Helianthus" from TS-Optics Set Dauerpräparate Botanik 25St. Objective: EC Plan-Neofluar 10x/0.3

Your Insight into the Technology Behind It

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service



Boost the speed of SIM imaging even further

Elyra 7 already provides you with fast imaging speeds. But you can further increase the temporal resolution and productivity for 2D and 3D imaging by using the speed enhancement modes. The Burst mode and the Leap mode are compatible with Lattice SIM as well as SIM Apotome acquisition. Combined with SIM² image reconstruction, they enable you to capture highly dynamic processes at exceptional resolution in all three dimensions.

2D Burst mode:

Get full temporal information

Burst mode processing uses the rolling window approach to let you observe processes in your living samples at up to 255 fps. Since Burst mode is a post-acquisition step, you have the flexibility to use it with previously acquired data sets. You decide how much temporal resolution is required for your data analysis.





U2OS cell expressing an Rab5-mEmerald (green) and tdTomato tagged Golgi associated transport marker (orange). Simultaneous dual-color acquisition. Objective: Plan-Apochromat 63×/1.4 Oil

3D Leap mode:

Digital sectioning at a new level

For demanding fast imaging in 3D, the Leap mode acquisition enables you to reduce your imaging time and lower the light exposure on your sample. This works by imaging only every third plane, for three-times higher volume imaging speed and three-times fewer light exposures. ZEN reconstructs the entire volume using a pixel reassignment approach.



Cos-7 cells stained with anti-beta-tubulin Alexa Fluor 488. Images show XY, XZ and YZ view of the cropped volume image for Nyquist sampled (left) and digitally sectioned (right) images of the same sample area. Objective: Plan-Apochromat 63×/1.4 Oil

Expand Your Possibilities




- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

SMLM:

3D imaging at molecular resolution

In single-molecule localization microscopy (SMLM), fluorescent molecules are sparsely activated so that only one out of many will be in its on-state within a single point spread function (PSF). This lets you determine its center of mass with a localization precision that far exceeds the extension of the PSF. Once recorded, the molecule is turned to its off-state and the cycle of activation/deactivation is repeated until all molecules are captured. The localizations are plotted in a new image to create the super-resolution image. With Elyra 7 you can use SMLM techniques such as PALM, dSTORM and PAINT to achieve lateral localization precision of 10–20 nm. The ZEN software will seamlessly perform the image reconstruction of your data.

In addition, Elyra 7 provides you with 3D SMLM mode based on PRILM technology. The PSF is reshaped for encoding the Z position so while acquiring only one plane, you get volume information of 1.4 μm depth at 20–40 nm axial resolution. Thus, you can acquire 3D data with consistent molecular precision.



3D PAINT image of mitochondrial membranes in BSC1 (kidney epithelial cells). The outer membrane protein Tomm20 was labeled using Ultravue – I2-650 imaging strand. Reshaped PSF encoding for Z information was used to create a 1.4 μm deep 3D PAINT image. Upper figures show the widefield microscopy image (left) of the area as well as the 3D image of the modified PSF (right) used for 3D SMLM. Objective: alpha Plan-Apochromat 63 \times /1.46 Oil

Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service


Simultaneous two-color imaging

Investigation of living samples very often focuses on interactions of different proteins or organelles. Simultaneous imaging of the involved structures is key to proper understanding of these highly dynamic processes. Equip your ZEISS Elyra 7 with a Duolink adapter to operate two Hamamatsu ORCA-Fusion sCMOS cameras in parallel and use all the advantages that a widefield-based technology can offer:

- True simultaneous two-color imaging within your entire field of view – without any delays
- Acquisition of a super-resolved real-time snapshot of an entire living cell by picking a low exposure time
- Increased productivity of fixed cell experiments by doubling the information obtained at the same time
- Imaging of any possible color combination, with minimal signal crosstalk as enabled by the integrated multi-bandpass emission filters
- Acquisition of 4-color images without mechanical filter change – making your multi-color experiments even faster
- Multi-color SMLM experiments

Hamamatsu ORCA-Fusion BT

This camera features a scientific CMOS (sCMOS) with a back-thinned sensor enabling peak quantum efficiency of ~95 %. With its ultra-low, uniform readout noise and CoaXPress interface for high acquisition speeds and exposure times down to 1 ms, it yields unparallel digital imaging results.



Cos-7 cell expressing the endoplasmic reticulum marker Calreticulin-tdTomato (magenta) and mitochondrial marker Tomm20-mEmerald (green) was simultaneously imaged for two colors. The movie shows high dynamic interactions of the ER and mitochondria. Objective: Plan-Apochromat 63x/1.4 Oil



Elyra 7 Duolink sCMOS camera adapter for simultaneous two-color acquisition with integrated multi-bandpass emission filter cubes for efficient image acquisition


Expand Your Possibilities

- › In Brief
- › **The Advantages**
- › The Applications
- › The System
- › Technology and Details
- › Service

Translate your images into quantitative data

ZEISS ZEN: A journey through different scales


Biological samples often contain different types of information at different length scales. Collecting low to high resolution data in the same sample not only makes you more productive, but also allows you to interconnect your findings and create more accurate biological models based on your experimental findings. With AI Sample Finder, you automatically detect your whole sample even before starting your experiment – ensuring that you won't miss any relevant areas. The ZEN Connect toolkit enables you to combine different experiments recorded with various acquisition modes or systems – placing your experiments into the spatial context of the whole sample.



Driven by ZEN (blue edition), ZEISS Elyra 7 is now fully integrated into the ZEN software ecosystem, providing workflow-oriented solutions for image acquisition, AI-based data processing, and advanced correlative analyses.

ZEISS arivis Pro: Advanced image processing and 3D reconstruction

Use the efficient ZEISS arivis Pro software for visualization and quantification of large 3D and 4D data sets. ZEISS arivis Pro not only renders volume images of almost unlimited size, but also provides advanced image processing tools such as volume fusion, channel shift, conventional and machine learning based segmentation, 3D tracking, and neuron tracing. Visualize your quantitative results within ZEISS arivis Pro or export all data for further analysis. The modular structure of ZEISS arivis Pro flexibly adjusts to your needs for advanced image processing and analysis.



Murine brain expressing the neuronal marker Thy1-eGFP, imaged in SIM Apotome and Lattice SIM modes over a Z stack range of 170 μm . The ZEN Connect project combines data sets recorded with 10 \times SIM Apotome, 25 \times SIM Apotome, 40 \times SIM Apotome and 63 \times Lattice SIM. The volume rendering shows a subset of the 63 \times Lattice SIM data set. Objective: Plan-Apochromat 63 \times /1.4 Oil. Sample courtesy of Herms Lab (MCN, University of Munich, Germany).

ZEISS Elyra 7 at Work



- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Lattice SIM² – Simply Image More



Study of the components of the cytoskeleton is a prominent research field in biology. Due to the fine structures of these components, for example the actin network or microtubular filaments, imaging far below 100 nm is often performed with super-resolution techniques.

Lattice SIM² now allows you to gain much more structural information from your samples compared to conventional SIM techniques. It not only operates with a resolution of down to 60 nm but also provides markedly improved sectioning quality in your images. This novel, robust image reconstruction method efficiently separates signal and background – all without the need for tailor-made staining protocols or expert knowledge of complex microscopy techniques.

You can take advantage of the easy-to-use Lattice SIM² technology to unveil complex structural information and simply see more.



The Lattice SIM² images of Cos-7 cells labeled with phalloidin Alexa 488 were acquired with an alpha Plan-Apochromat 100×/1.57 oil objective. Maximum intensity projection of Z stack is shown.



The Lattice SIM² images of Cos-7 cells labeled via immunofluorescence with anti-alpha-Tubulin Alexa 488 are shown as color-coded projection. Data were acquired with an alpha Plan-Apochromat 100×/1.57 oil (left) and Plan-Apochromat 63×/1.4 oil (right) objectives. The images demonstrate the excellent sectioning capabilities of SIM² image reconstruction algorithm.


ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Straightforward multi-color super-resolution imaging

Studying multiprotein complexes requires super-resolution imaging with multiple colors, which is not often attainable with conventional techniques. Lattice SIM² enables you to perform multi-color imaging at resolution down to 60 nm for conventionally stained samples.

The synaptonemal complex is a well-known structure in the nucleus of meiotic cells and consists of two lateral elements, which are connected to a central element by transverse filaments. Due to its small size, three-color imaging of the synaptonemal complex has previously been possible only using complex methods with elaborate sample preparation, such as super-resolution imaging of three-fold expanded samples using the expansion microscopy technique. Lattice SIM² resolves the two strands of SYCP3 (lateral elements) as well as SYCP1-C (C-terminus of transverse filaments) without special sample treatment or staining for distances well below 100 nm. More importantly, the three-color image provides structural information for the distances between the proteins SYCP3 and SYCP1. Even within the SYCP1 protein, the differently labeled N- and C-Terminus can be clearly separated with less than 50 nm resolution between the two labels.



Architecture of threefold labeled synaptonemal complexes from mouse testis visualized via immunolabeling of SYCP3 with SeTau647, SYCP1-C with Alexa 488 and SYCP1-N with Alexa 568 and Lattice SIM² mode. Sample courtesy: Marie-Christin Spindler, AG Prof Ricardo Benavente, Biocenter of the University of Würzburg. Objective: Plan-Apochromat 63×/1.4 Oil


ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Observe life's finest details

ZEISS Elyra 7 was designed for the understanding of biological processes that require the observation of living samples in low light and high spatiotemporal resolution. Observe cellular, sub-cellular, and sub-organelle structures in living specimens in 2D and 3D over time. Whether you are interested in cytoskeletal dynamics, mitochondrial fusion and fission, or budding of the endoplasmic reticulum, it gives you both live cell compatibility and super-resolution.


To image the fine structures of the cytoskeleton components, e.g., actin and ER networks or microtubular filaments, super-resolution imaging below 100 nm is required. Mitochondria are powerhouses of our cells, generating energy in the form of ATP to sustain the cell. These highly dynamic organelles constantly undergo fusion and fission events to ensure proper distribution of ATP across the cell. They are known to interact with the cytoskeleton and other subcellular compartments, e.g., microtubules, which they ride on to get to their destinations, or the ER, which wraps around mitochondria to initially constrict their diameter before fission events.



Actin dynamics in a U2OS cell expressing LifeAct-GFP were imaged with the Lattice SIM 3D Leap mode and reduced phases. The movie shows a color-coded depth projection of the volume stack. Objective: Plan-Apochromat 63×/1.4 Oil



U2OS cell expressing Tomm20-mEmerald. Image shows a color-coded projection of the Lattice SIM² volume data set. Objective: Plan-Apochromat 63×/1.4 Oil




Simultaneous imaging of the endoplasmic reticulum (Calreticulin-tdTomato, magenta) and microtubules (EMTB-3xGFP, green) in a Cos-7 cell reveals highly dynamic interaction of these organelles. Objective: Plan-Apochromat 63×/1.4 Oil

ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Excellent sectioning at incredible speed

SIM² Apotome is your flexible live cell imaging method for experiments that do not require the highest spatial resolution but rely instead on excellent sectioning quality. SIM² Apotome is superior to conventional confocal microscopy in terms of lateral and axial resolution as well as volume acquisition speed while it is also very gentle to your sample. Here, LLC PK1 cells expressing H2B-mCherry and α -Tubulin mEmerald-GFP were continuously observed with a 25 \times /0.8 water immersion objective over 35 min, while undergoing mitosis.



SIM² Apotome time lapse data of LLC PK1 cells expressing H2B-mCherry (magenta) and α -Tubulin mEmerald-GFP (green). Data shown as maximum intensity projection of 12 planes over 3.7 μ m depth. Objective: LD LCI Plan-Apochromat 25 \times /0.8 Imm Corr

SIM² Apotome is compatible with objectives of different magnifications (10 \times , 20 \times , 25 \times and 40 \times). The high NA (1.4) 40 \times magnification images almost reach the resolution and sectioning capabilities of a conventional SIM microscope as demonstrated for Cos-7 cell expressing Calreticulin-tdTomato, while multiplying acquisition speed.




SIM² Apotome time lapse data of Cos-7 cells expressing the endoplasmic reticulum marker Calreticulin-tdTomato. Data shown as maximum intensity projection of 12 planes over 1.4 μ m depth. Objective: Plan-Apochromat 40 \times /1.4 Oil

ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Super-resolution imaging at up to 255 fps


The diffusive and especially the ballistic movement of small vesicles in cells can be captured only when super-resolution and high-dynamic imaging are possible at the same time. With the Burst processing of 2D time lapse data, Elyra 7 is able to generate super-resolution images at 255 Hz in a large field of view and even acquire two colors simultaneously in both Lattice SIM and SIM Apotome acquisition modes.



U2OS cell expressing Rab5-mEmerald (green) and tdTomato tagged Golgi associated transport marker (magenta). Simultaneous dual-color acquisition with an exposure time of 1.5ms/phase for a FOV of 1024×1024 pixel (64 μm×64 μm). Objective: Plan-Apochromat 63×/1.4 Oil

Digital sectioning for 3D imaging three times faster

Elyra 7 Leap mode accelerates the volume imaging speed three times and at the same time decreases the light dosage on your sample. While still capturing all the finest details, the entire volume (18 planes) of the U2OS cell expressing Calreticulin-tdTomato was imaged at 38 volumes/min speed in Lattice SIM acquisition mode. For SIM Apotome acquisition mode, you can expect up to three times higher volume imaging speed.




U2OS cell expressing calreticulin-tdTomato to visualize the endoplasmic reticulum. The time series shows a maximum intensity projection of the volume data set. Objective: Plan-Apochromat 63×/1.4 Oil

ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Resolve the details hiding in the depth

Despite being a structured illumination-based microscope, Elyra 7 Lattice SIM² as well as SIM² Apotome also provide you with super-resolution and high-quality sectioning in thick or scattering samples.



Murine brain expressing the neuronal marker Thy1-eGFP was imaged in Lattice SIM mode over a Z stack range of 75 μm. The image shows the color-coded projection of the volume data. Objective: Plan-Apochromat 63×/1.4 Oil. Sample courtesy of Herms Lab (MCN, University of Munich, Germany)

The combination of robust illumination patterns and excellent image reconstruction technology enabled us to image throughout an entire murine brain section of ~80 μm thickness expressing the neuronal marker Thy1-eGFP.



Zebrafish embryo expressing a vascular marker fli1-EGFP was imaged in SIM Apotome mode over a Z stack range of 100 μm. The SIM² processed image shows the color-coded projection of the volume data. Objective: LD LCI Plan-Apochromat 25×/0.8 Imm Corr. Sample courtesy of Haass Lab (MCN, University of Munich, Germany)

ZEISS Elyra 7 at Work


- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service

Correlative microscopy within the same system


Correlative microscopy, where the same area of the sample is imaged using different techniques, has become an important tool for biological research. Imaging with Elyra 7 can be combined seamlessly with LSM 980 Airyscan or ZEISS electron microscopy solutions. Nevertheless, Elyra 7 itself contains three different imaging modalities – SIM Apotome, Lattice SIM and SMLM – providing the possibility to combine them when needed for sample scales varying by orders of magnitude.

Here, a detailed SIM² Apotome overview volume image of a drosophila larvae was combined with the super-resolution Lattice SIM² volume image of a region of interest, putting the super-resolved area into full context.


Another advantageous combination is the pre-imaging of SMLM samples with Lattice SIM² to easily identify interesting sample areas at resolutions of 60–100 nm, then perform the more time-consuming localization microscopy on suitable cells.



SIM² Apotome and Lattice SIM² images of *D. melanogaster* larva stained with HLH-54F-GFP (cyan) and Anti-Cut-Cy3 (magenta). Images show maximum intensity projections of 3D data. Sample courtesy of R. Palmer and G. Wolfstetter (University of Gothenburg)



Lattice SIM² and dSTORM images of bacteria stained with a membrane marker coupled to Alexa Fluor 647. Sample courtesy of J. Nabarro, C. Baumann, G. Calder & P. O'Toole (Department of Biology & Bioscience Technology Facility, University of York, UK)



ZEISS Elyra 7 at Work

- › In Brief
- › The Advantages
- › **The Applications**
- › The System
- › Technology and Details
- › Service


Single-molecule localization microscopy (SMLM)

SMLM encompasses techniques such as PALM, dSTORM, and PAINT. With high power lasers across the visible spectrum and dual camera detection, Elyra 7 allows researchers to gain access to a broad range of dyes, markers and fluorophores in almost any possible combination.

Elyra 7 enables precise quantification over a large field of view and an unprecedented Z-capture range. You now can acquire 3D data from a whole cell with molecular precision.

Resolve molecular structures


SMLM allows you to map precise locations of individual proteins.



SMLM: Eightfold symmetry of the nuclear pore complex in A6 cells. Gp210 was labeled with Alexa Fluor 647. Widefield image (top left), SMLM image (top right) and zoomed in region (bottom).

Determine the relationships between molecules


Detect two channels with molecular precision.



Dual-color 2D STORM of Cos-7 cells stained for microtubules (anti-tubulin-Alexa Fluor 647, magenta) and mitochondria (anti-TOMM20-CF568, yellow).

Capture information in three dimensions

Untangle molecular relationships in Z with confidence.



SMLM: With Elyra 7 you can image a Z-depth of 1.4 μm in a single acquisition. 3D SMLM image of Alexa 647 α-tubulin color-coded for depth.

Sample courtesy of Michael W. Davidson, Florida State University, USA.

The Lattice SIM product family

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service

Address all your super-resolution needs across scales

The ZEISS Lattice SIM product family gives you full access to super-resolution imaging for all research areas, from fast optical sectioning to the detection of highly dynamic processes and quantification at the molecular level.



ZEISS Lattice SIM 3

Reveal cellular behavior and inter-cellular dynamics

Lattice SIM 3 is specifically designed to meet the imaging requirements of multicellular organisms and tissue sections. This system exploits the full potential of the SIM Apotome technology: fast optical sectioning at superior quality, large fields of view with access to smaller regions of interest, near-isotropic resolution, and the gentlest super-resolution imaging possible.



ZEISS Lattice SIM 5

Reveal the vibrant sub-organelle network of life

ZEISS Lattice SIM 5 has been optimized for single cell imaging as well as capturing subcellular structures and their dynamics. Powered by the Lattice SIM technology and the SIM² image reconstruction algorithm, ZEISS Lattice SIM 5 provides you with outstanding super-resolution capabilities down to 60 nm in both living and fixed cells.




ZEISS Elyra 7 with Lattice SIM

Reveal life across scales – down to molecular details

ZEISS Elyra 7 includes several microscopy techniques: Lattice SIM², SIM² Apotome, SMLM and TIRF. You can combine these techniques to multiply the insights from your specimen and to correlate the acquired data. With its focus on single molecule localization microscopy, ZEISS Elyra 7 gives you resolution excellence down to the molecular level.

Your Flexible Choice of Components

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



1 Microscope

- Axio Observer 7 (inverse stand)
- Stage top incubation
- Motorized XY stepper scanning stage
- Z-Piezo stage insert
- 2 camera ports or one camera port with Duolink

2 Objectives

- C-Apochromat 63×/1.2 Water (DIC)*
- Plan-Apochromat 63×/1.4 Oil (DIC)*
- alpha Plan-Apochromat 100×/1.46 Oil (DIC)*
- alpha Plan-Apochromat 100×/1.57 Oil HI Corr (DIC)*
- alpha Plan-Apochromat 63×/1.46 Oil
- C-Apochromat 40×/1.2 W
- Plan-Apochromat 40×/1.4 Oil (DIC)*
- LD LCI Plan-Apochromat 25×/0.8 Imm Corr
- Plan-Apochromat 20×/0.8 Air
- EC Plan-Neofluar 10×/0.3 Air

3 ZEISS Elyra 7 Illumination and Detection

- Fiber coupled solid state or diode pumped solid state lasers
- Available lines:
 - 405 nm diode (50 mW),
 - 488 nm OPSL (100 or 500 mW),
 - 561 nm OPSL (100 or 500 mW),
 - 642 nm diode (100 mW),
 - 640 nm OPSL (500 mW)
- Lasers shared between Lattice SIM and SMLM
- Hamamatsu ORCA-Fusion BT sCMOS camera


4 Software

- ZEN (blue edition)
- SIM toolkit
- SMLM (PALM/dSTORM) module

* DIC indicates type of objective not imaging modality


System Overview

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Incubation Setup

- › In Brief
- › The Advantages
- › The Applications
- › **The System**
- › Technology and Details
- › Service



Technical Specifications

- » In Brief
- » The Advantages
- » The Applications
- » The System
- » **Technology and Details**
- » Service

Microscope

Stand	Axio Observer 7, motorized inverted microscope for superresolution microscopy
Z-drive	DC servomotor, opto-electronically coded; smallest Z step 25 nm
XY Piezo scanning stage	motorized; range 130 mm × 100 mm; max sped 100 mm/s; resolution 0.2 µm; reproducibility: ± 1 µm; absolute accuracy ± 5 µm; suitable for mounting frames K 160 × 110 mm and Z-Piezo Stage insert
Z-Piezo stage insert	for XY scanning stage, max travel range 500 µm; smallest Z step size 5 nm, sample holders available for standard 3"×1" slides LabTek chambers, multiwell plates and 36 mm glass-bottom dishes; level-adjustable and universal stage insert available for standard slides, glass-bottom dishes and LabTek™ chambers.

Optical Filters for Lattice SIM and SMLM

Filter sets reflector turret	Four exchangeable filter sets available for multi-channel Lattice SIM and SMLM; each filter set with four precisely mounted ACR-coded ⁽¹⁾ filter modules for superresolution microscopy on a motorized six-position turret; two positions in each turret compatible with standard Push & Click filter modules, e.g. for visual sample observation.
Dual filter sets for Duolink optimized for dual color and double dual color applications	Filter sets are optimized for dual camera applications, maximum sensitivity, minimal cross-talk and reduced autofluorescence.
Filter slider	Manual filter slider with two positions (for emission filters or a Bertrand lens); fits into camera adapter of the microscope's side port; emission filters exchangeable for customizing detection conditions.

Lasers

Laser module for Elyra 7	Laser coupling with polarization-maintaining single mode fiber (no adjustment of laser coupling by users required).
Laser lines	405 nm (50 mW), 488 nm (100 mW or 500 mW), 561 nm (100 mW or 500 mW), 642 nm (100 mW) or 640 nm (500 mW); 405 laser can be attenuated by up to 100000 fold (used for activation and back-pumping); high power lasers (500 mW) can be 10 fold attenuated (488, 561, 642)

Cameras

sCMOS	Hamamatsu ORCA-Fusion BT; sensor pixel count: 2304 × 2304, effective: 1288 × 1288; pixel size 6.5 µm × 6.5 µm; QE: up to 95 % @ 540 nm); water cooling (regulated sensor temperature −8 °C); dynamic range: 16 bit; binning: 1 × 1, 2 × 2, 4 × 4; Frame rates: 89 fps (fast) @ full frame in Laser WF modality Liquid cooling system for sCMOS camera (Hamamatsu ORCA-Fusion BT)
-------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Elyra 7 for SMLM

Illumination module	<p>Type L and LS: Fully motorized Epifluorescence (EPI), high inclined and laminated optical sheet (HILO) and total internal reflection illumination (TIRF); simultaneous TIRF illumination with VIS and 405 nm laser lines; individual triggering of lasers for synchronizing dye activation and illumination to camera read-out and transfer times; motorized TIRF angle adjustment; motorized TIRF field adjustment with three field size options</p> <p>Type LS: Five different grating frequencies for Lattice SIM for optimal matching of illumination pattern to laser wavelength and objective lens; Motorized exchange of gratings in multi-color Lattice SIM; one grating for SIM Apotome; fast piezo actuated phase stepping of gratings</p>
3D-PALM module	Double phase ramp in pupil plane of back aperture of objective providing for phase ramp imaging localization microscopy (PRILM); z capture range typically 1.4 μm
Camera	Up to two Hamamatsu ORCA-Fusion BT sCMOS cameras mounted on right side port
Objective lenses (SMLM)	alpha Plan-Apochromat 100 \times /1.46 Oil DIC*, alpha Plan-Apochromat 100 \times /1.57 Oil-HI DIC* Corr (2D-PALM), alpha Plan-Apochromat 63 \times /1.46 Oil, Plan-Apochromat 63 \times /1.4 Oil DIC*, C-Apochromat 63 \times /1.2 W Corr DIC* (3D-PALM) ACR ⁽¹⁾ coding (optional); Objectives with NA \geq 1.46 suitable for TIRF and HILO illumination)
Imaging modes	Widefield (WF) mode (sample illumination with thermal light source or LED), Laser WF mode (sample illumination with laser), SMLM mode for single-molecule localization microscopy
Field of view (SMLM)	Maximal field of view 50.9 \times 50.9 μm^2 (with alpha Plan-Apochromat 100 \times /1.46 Oil DIC*, 1.6 \times tube lens, full chip recording); 128.22 \times 128.22 μm^2 (with Plan-Apochromat 63 \times /1.4 Oil DIC*, 1.0 \times tube lens, full chip recording); HP field 2 \times smaller, uHP field 2 \times $\sqrt{2}$ smaller than TIRF field
Localization precision (SMLM)	Typically 10 nm – 20 nm lateral, 20 nm – 40 nm axial, given sufficient signal-to-noise & density
Multi-color imaging (SMLM)	Detection of up to two different fluorescent labels (simultaneous with Duolink or quasi simultaneously by fast sequential laser switching)
Acquisition speed (SMLM)	sCMOS (dSTORM) and widefield mode > 100 frames per second full chip; 200 frames per second (512 \times 512 pixels)
Data recording and analysis (SMLM)	<p>Full software control of SMLM imaging; Definite Focus z-drift control</p> <p>Online SMLM processing for simultaneous data acquisition and analysis; manual editing of parameter settings for optimal results in SMLM with different fluorophores; feature-rich rendering of SMLM localization tables; export and import of localization tables for custom filtering; correction algorithms for lateral and axial drift; chromatic aberration correction (based on fiducial markers or prominent structures)</p> <p>Multi-emitter fitting algorithms allow to analyze overlapping signals with high precision. Up to 10 times higher labeling densities are possible speeding up acquisitions by the same factor.</p>

* DIC indicates type of objective not imaging modality

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Elyra 7 for Lattice SIM and SIM Apotome mode

Illumination module	Fully motorized Lattice SIM imaging; five different grating frequencies for Lattice SIM for optimal matching of illumination pattern to laser wavelength and objective lens; motorized exchange of gratings in multi-color Lattice SIM; one grating for SIM Apotome fast piezo actuated phase stepping of gratings.
Camera	Up to two Hamamatsu ORCA-Fusion BT sCMOS cameras mounted on right side port
Imaging Modes	Widefield modes for illumination with thermal light source or LED and lasers, Lattice SIM using two dimensional grid SIM mode (two- and three-dimensional Lattice SIM), SIM Apotome mode using one dimensional grid for z-sectioning
Objective lenses (Lattice SIM)	Plan-Apochromat 63×/1.40 Oil DIC*, C-Apochromat 63×/1.20 W Corr, alpha Plan-Apochromat 63×/1.46 Oil, ACR ⁽¹⁾ coding (optional), alpha Plan-Apochromat 100×/1.57 Oil-HI DIC* Corr
Objective lenses (SIM Apotome mode)	Plan-Apochromat 40×/1.4 Oil; C-Apochromat 40×/1.2 W; EC Plan-Neofluar 10×/0.3 Air; Plan-Apochromat 20×/0.8 Air; LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*
Resolution (Lattice SIM)	Lateral resolution (XY): down to 120 nm, axial resolution (Z): down to 300 nm (typical experimental FWHM values with objective lens Plan-Apochromat 63×/1.40 Oil DIC*, subresolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Resolution (Lattice SIM ²)	Lateral resolution (XY): down to 60 nm, axial resolution (Z): down to 200 nm (typical experimental FWHM values with objective lens Plan-Apochromat 63×/1.40 Oil DIC*, subresolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Resolution (SIM ² Apotome)	Lateral resolution (XY) down to 140 nm for 40× 1.4 oil objective; lateral resolution (XY) down to 265 nm for 25× 0.8 Imm Corr objective; (typical FWHM measurements on sub-resolution beads of 100 nm diameter and excitation at 488 nm; resolution is sample and SNR dependent)
Multi-color (Lattice SIM and SIM Apotome mode)	Detection of up to four different fluorescent labels (sequential detection) and simultaneous dual color detection with DuoLink
Max. Field of view (Lattice SIM)	80.14×80.14 μm ² , full-frame recording (1280×1280 effective px) with Plan-Apochromat 63×/1.40 Oil DIC*
Max. Field of view (SIM Apotome mode)	126×126 μm ² , full frame recording (1280×1280 effective px) with Plan-Apochromat 40×/1.40 Oil 203×203 μm ² , full frame recording with LD LCI Plan-Apochromat 25×/0.8 Imm Corr DIC*; 252×252 μm ² , full frame recording with Plan-Apochromat 20×/0.8 Air; 505×505 μm ² , full frame recording with EC Plan-Neofluar 10×/0.3 Air
Acquisition speed (Lattice SIM)	19 SIM image frames per second at 512×512 px resolution and 1 ms exposure time (13 phase images per one SIM image) 28 SIM image frames per second at 512×512 px resolution and 1 ms exposure time (9 phase images per one SIM image)
Acquisition speed (SIM Apotome mode)	51 sectioned frames per second at 512×512 px resolution and 1 ms exposure time (camera limited) (5 phase images per one sectioned image); 85 sectioned frames per second at 512×512 px resolution and 1 ms exposure time (camera limited) (3 phase images per one sectioned image);
Leap mode and Burst mode	Leap and Burst modes are combinable with both the Lattice SIM and SIM Apotome. Leap mode increases the frame rate by a factor of 3 for 3D image acquisition. Max. 255 image frames per second at 512×512 px resolution and 1 ms exposure time are available for 2D data after Burst processing.
Data recording and analysis (Lattice SIM and SIM Apotome mode)	Full software control of Lattice SIM imaging; Multi-tracking (sequential multi-channel data acquisition with freely configurable change of gratings (Lattice SIM), or one common grating (SIM Apotome mode), filters and excitation lasers between tracks); Simultaneous dual color imaging with one grating; Lattice SIM and SIM Apotome mode imaging in user-defined sub-array regions (ROI imaging); Leap mode for 3 times faster imaging with excellent sectioning; Extension of imaged area possible with tile scanning and stitching. Burst mode processing for 2D time series data sets for Lattice SIM and Aptotome mode to increase effective frame rates by a factor of 15 and 5, respectively.

* DIC indicates type of objective not imaging modality

⁽¹⁾ ACR (Automatic Component Recognition); Elyra 7 systems and ZEN imaging software automatically recognize ACR-coded components.

Technical Specifications

- › In Brief
- › The Advantages
- › The Applications
- › The System
- › **Technology and Details**
- › Service

Elyra 7 for combined Lattice SIM and SMLM

System information	All imaging modes combined in one system
Illumination module	Sample illumination in all widefield and superresolution modes by a single, highly integrated illumination module (with same set of lasers and a single Elyra laser module).
Camera	Up to two Hamamatsu ORCA-Fusion BT sCMOS cameras mounted on right side port

Software

Standard	ZEN imaging software (64-bit); operating system: Microsoft Windows 10 Full software control of image data recording in all imaging modes (including widefield, superresolution); Software-controlled switching between imaging modes. Full software control of data recording (multi-channel imaging, time series, z-stack) Saving and restoring of user-specific configurations for data recording.
Software packages	Required: ZEN package; ZEN module Lattice SIM; ZEN toolkit Advanced Acquisition; ZEN toolkit 3D Optional: ZEN toolkit 2D; ZEN toolkit Deconvolution; ZEN toolkit Connect; ZEN toolkit AI; ZEN toolkit Developer; Vision package

Accessories

Definite Focus	Holding focus to compensate axial drift, typical z-position accuracy: 30 nm; Specified limits of Definite Focus 3: $0.2 \times \text{DOF}$ (Depth of field: $\text{DOF} \approx \lambda / \text{NA}^2$)
Incubation	Stage-top incubation possible with safety lock
Duolink for attachment of two cameras of the same type	Allows attachment of two cameras of the same type to the microscope.
Storage PC with 55 TByte storage capacity	Direct streaming of data and parallel processing while streaming of data possible



Elyra 7 with Lattice SIM meets the requirements according to IEC 60825-1:2014 and is a laser class 3B device. Interlocks on customer interfaces prevent access to the laser radiation.



По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: zsf@nt-rt.ru || сайт: <https://zeiss.nt-rt.ru/>